

24hrs: 0086-18021231901
www.saikangmedical.com

SAIKANG (Mid East Branch)
Saikang Medical Supplies (Mid East) LLC
Muscat
PO Box 53200, Abu Dhais, UAE
Tel: +971 585710652

SAIKANG (Russia Branch)
LLC "SAIKANG"
454004, Russian Federation, Chelyabinsk,
Universitetskaya Naberezhnaya
TEL: +7 (919) 300-00-06

SAIKANG (South America Branch)
Saikang Medical Technology, Inc.
W3400 Corrientes, Argentina
TEL: +55 (11) 2758-0438
Rua Universal, Hollywood,
São Bernardo do Campo, Brazil
TEL: +55 (11) 2758-0438

SAIKANG (East Africa Branch)
Saikang Medical Technology Ltd.
Mezzanine Floor, Taj Tower,
Upper Hill, Nairobi, Kenya
TEL: +254714084806

SAIKANG (South East Asia Branch)
Saikang Medical Technology Ltd.
Sala Thammasop Subdistrict, Thuan
Wattana District
Bangkok, Thailand
TEL: +66-642741046

SAIKANG (China Headquarter)
No.25 Leihong Road, Zhonggang City
Jiangsu Province, China
TEL: +86-512-58525850
E-mail: inquiry@saikang.com

SAIKANG (UK Branch)
Uk Healthward International Limited
DEPT 111,
196 High Road Wood Green,
London, N22 8HH
United Kingdom
TEL: +31 (0) 202111106



Factory Tour

Website



SAIKANG
Leading Brand of China Medical Furniture

www.saikangmedical.com
E22-1

SAIKANG®

LEADING BRAND OF CHINA
MEDICAL FURNITURE

JIANGSU SAIKANG MEDICAL EQUIPMENT CO., LTD.
www.saikangmedical.com

Stock Code: 870098



SAIKANG[®]
LEADING BRAND OF CHINA
MEDICAL FURNITURE



LEADING BRAND OF CHINA
MEDICAL FURNITURE

RACE AGAINST TIME
HEALTHIER YOUR LIFE



CONTENTS

01 Strength of Saikang	09 Hospital Bed	63 Mattress
03 Technology of Saikang	29 Homecare Bed	65 Bedside Table
05 Innovation of Saikang	33 Orthopedic Bed	66 Overbed Table
07 Quality of Saikang	35 Children Bed	67 Medical Chair
	39 Baby Crib	76 Stretcher
	40 Bathing Bed	78 Stretcher Trolley
	41 Examination Table	79 Patient Trolley
	44 Gynecological Table	85 Medical Trolley
	47 Obstetric Bed	91 Surgical Instrument
	51 Operating Table	92 Cabinet/Shelf

STRENGTH OF SAIKANG

20

20 years focusing on design, production, service and technical consultation of medical care equipment.

150

Located in Zhonggong City, 150km nearby Shanghai Economic Center.

10000000

Saikang delivered 1,000,000 sets of products globally. Products all over Asia, Middle East, North America and Africa.

400

There are more than 400 employees worldwide, more than 50 qualified engineers, and more than 100 trained technicians.

189

Equipped with 189 sets of various automatic production facilities, Saikang follow strict standard control system on every production process, enable the traceable of each part.

50000

The finished phase I production area covers more than 30,000 square meters, and phase II will be more than 20,000 square meters of advanced production area.

18

Saikang offer comprehensive medical solutions with 18 product categories.

562

Safe, efficient and easy-to-use medical solution including emergency nursing and ICU wards of more than 500 series of products.

29

Win CE certification, ISO9001/13485 quality system certification, U.S. FDA and other 29 certifications.

131

Ownership of 131 patents, continuously leading new standards in defining product quality and service efficiency.

149

Since the establishment of the company, products and solutions have been used in more than 149 countries.

1000

Saikang products are committed to creating a modern nursing environment, accumulatively serving more than 1,000 medical and health care organizations.



LEADING BRAND OF MEDICAL FURNITURE

Patients in medical rehabilitation often have to face huge challenges. Many people who suffering from diseases or accidents are faced with retraining their minds and bodies and must adapt to the new environment. In this case, the positive atmosphere would offer decisive help. Hospital beds and furniture should provide the best support for patients and nurses. Considering this case, Saikang devoted to provide medical care solutions and perfect choices for medical institutions and rehabilitation centers, to help patients achieve rapid or long-term recovery.

Saikang was founded in 2002, now is a leading global provider of medical furniture and solutions. Head office locates in Suzhou, and branches in U.S., U.K., Dubai, Russia and Kenya, having 400 employees worldwide. Saikang has formed a large global network of research and development, marketing and service. Currently, Saikang products and medical solutions have been used in more than 108 countries and regions in the world.

TECHNOLOGY OF SAIKANG

SAIKANG insists on providing high quality products. All products have passed the "Made in China" products in world's first-class. SAIKANG has ISO9001/14001 quality system certification, conforms to EU safety standards and has passed CE certification, ensuring the market access worldwide.

Alibaba.com Assessed Supplier

Jiangsu Saikang Medical Equipment Co., Ltd.
江苏赛康医疗设备有限公司

中国制造网 认证供应商
AUDITED SUPPLIER
Presented to
Jiangsu Saikang Medical Equipment Co., Ltd.
江苏赛康医疗设备有限公司

SGS
Made in China
中国制造商





INNOVATION OF SAIKANG

Design Innovation



SAIKANG medical designers have always maintain close communication with medical service providers, to help them providing better solutions to their patients. With long-standing experience in product structure and functional experience, SAIKANG continue to forge ahead with innovation and play a vital role in helping medical service providers to improve the quality of care.

Research Innovation



SAIKANG medical owns independent research and development center. Since its inception, every year SAIKANG invested 10% of its annual sales in R&D, launched 10 new products, and applied at least 10 patented technologies on an average basis, which lead SAIKANG to continuously define the quality of medical products and new standards of efficiency.

Manufacturing Innovation



SAIKANG medical, as a pioneer in the field of medical furniture, always choose branded suppliers strictly to ensure the safety and reliability of raw material. SAIKANG adopt robot welding, automatic epoxy coating, efficient assembly line, dedicated operator, rigorous quality control, only to guarantee SAIKANG's commitment on quality.

Maintenance Innovation



SAIKANG medical have more than 50 engineers, and qualified service team that composed of professional subcontractors, trained technicians, and after-sale staff. They are experienced and have excellent skills, not only to respond each service request on time, but also precaution every detail in equipment operation, always focus on products safety.

Service Innovation



SAIKANG medical have thorough service architecture: 7days x 24hours online customer call center, complete spare parts warehouse, advanced CRM system, precautionary service, etc. No matter where customer is, SAIKANG service will be always close at hand.

QUALITY OF SAIKANG

Bed panels

A Harmonious Whole

Different panel design available, easy to remove and fix, without any tools. Fitted with buffers that protect both the walls and the bed itself when travelling.

Sleeping Surface


Unrivaled Sleeping Comfort

With 4 different sleeping surface to choose from: Profiled steel, thermoformed synthetic, shaped Plastic.

Protection Side Rails

Safety is Priceless

Atmosphere galleries, ABS collapsible side rails offer the patient protection on 2/3 of the bed. Thanks to the special design, the galleries constitute the ultimate protection for the patient.



Control Panel

Autonomous Control

Functional and easy to handle. The control offers all the functions which allow the patient and staff to choose according the positions.

Height Adjustment System

Perfect Stability

The height-adjustable support arms guarantee excellent bed stability. Thanks to the special design [patented] the patient can more easily get out of bed and receive care.

World Class Branded Parts

Perfect Partners

We support and produce our products in accordance to EN international standard, our outsourced parts are from world class tracked companies.

With a total area of more than 30,000m², SAIKANG can fulfill various requirements of production. Through the intelligent management system, each process is visualized and standardized. The manufacturing base has its strict quality management system, so every production procedure is enabled to follow their standards and every produced piece is traceable.



Technique. Be Unique.



Quality. Be Confident



Design. Be Harmonized



Oasis Series

Y8y

Electric ICU Bed
SKD-C

Technical Parameters

External Length	2130mm
External Width	1050mm
Bed Height Min.	360mm
Bed Height Max.	780mm
Mattress Platform	1925×880mm
Safe Working Load	250kg

Technical Configuration

✓ Linear Motor	4pcs
✓ 5" Double Side Castors, Central Locking System	4pcs
✓ ABS Oss Side Rail Extra Height	1set
✓ ABS Oss Bed Ends with Safe Lock	1set
✓ PP Bed Platform, Removable & Washable	9pcs
✓ Plastic Socket Inside PP Guardrail for Phone or Glasses	2pcs
✓ Patent Angle Indicator For Back-rest And Trendelenburg Position	4pcs
✓ Ø70mm Wall Bumper Wheel	4pcs
✓ IV Pole 4 Hooks	1set
✓ IV Pole Prevision	4pcs
✓ Drainage Hook	2pcs
✓ Central Locking Pedal	2pcs
✓ ABS Mattress Stopper	2pcs
✓ Hand Controller	1pc
✓ Manual CPR	
✓ Electric CPR	

Advantages

- ★ Double Auto-regression
- ★ Super-low Bed Position 350mm
- ★ PP Bed Platform, Removable & Washable
- ★ Inserted-type PP Head Foot Board, One Step Release With Safe Lock, Patent Design
- ★ PP Safety Guardrail With Two Steps Unlocking System, Patent Design
- ★ Back Rest 30degree Position, Reduce The Chance Of Lung Infection
- ★ Built-in LED Torch

Optional

<input type="checkbox"/> Weighing System
<input type="checkbox"/> Bed Extension 20cm
<input type="checkbox"/> Back-up Battery
<input type="checkbox"/> Nurse Control Panel
<input type="checkbox"/> Emergency Stop Button
<input type="checkbox"/> Cardiac Chair Position
<input type="checkbox"/> Directional Wheels Function
<input type="checkbox"/> 5th and 6th Extra Side Rail For Full Protection
<input type="checkbox"/> Monkey Pole
<input type="checkbox"/> Balkan frame
<input type="checkbox"/> Linen Holder
<input type="checkbox"/> One-button for patient mobilisation position
<input type="checkbox"/> One button nursing position
<input type="checkbox"/> Lowest position indicator



360° View



Oasis Series

Y8t

Electric ICU Bed
SKD-C

Technical Parameters

External Length	2140mm
External Width	1060mm
Bed Height Min.	480mm
Bed Height Max.	770mm
Mattress Platform	1925×900mm
Safe Working Load	250kg

Technical Configuration

✓ Linear Motor	4pcs
✓ 5" Double Side Castors, Central Locking System	4pcs
✓ ABS Oss Side Rail Extra Height	1set
✓ ABS Oss Bed Ends with Safe Lock	1set
✓ Stamping Forming Bed Platform	1pc
✓ Plastic Socket Inside PP Guardrail for Phone or Glasses	2pcs
✓ Patent Angle Indicator For Back-rest And Trendelenburg Position	4pcs
✓ Ø70mm Wall Bumper Wheel	4pcs
✓ IV Pole 4 Hooks	1set
✓ IV Pole Prevision	4pcs
✓ Drainage Hook	2pcs
✓ Central Locking Pedal	2pcs
✓ Emergency Stop Button	1pc
✓ Manual CPR	
✓ Electric CPR	

Optional

<input type="checkbox"/> Hand Controller
<input type="checkbox"/> Weighing System
<input type="checkbox"/> Nurse Control Panel
<input type="checkbox"/> Cardiac Chair Position
<input type="checkbox"/> Bed Extension 20cm
<input type="checkbox"/> Back-up Battery
<input type="checkbox"/> Directional Wheels Function
<input type="checkbox"/> ABS Mattress Stopper
<input type="checkbox"/> Monkey Pole
<input type="checkbox"/> Balkan Frame
<input type="checkbox"/> Linen Holder

Oasis Series

Y6y

Electric Bed
SKD-C

Technical Parameters

External Length	2130mm
External Width	1050mm
Bed Height Min.	360mm
Bed Height Max.	800mm
Mattress Platform	1925×880mm
Safe Working Load	250kg

Technical Configuration


✓ Linear Motor	3pcs
✓ 5" Double Side Castors, Central Locking System	4pcs
✓ ABS Oasis Side Rail Extra Height	1 set
✓ ABS Oasis Bed Ends with Safe Lock	1 set
✓ PP Bed Platform, Removable & Washable	9pcs
✓ Plastic Socket Inside PP Guardrail For Phone Or Glasses	2pcs
✓ Patent Angle Indicator For Back-rest And Trendelenburg Position	4pcs
✓ Ø70mm Wall Bumper Wheel	4pcs
✓ IV Pole 4 Hooks	1 set
✓ IV Pole Prevision	4pcs
✓ Drainage Hook	2pcs
✓ Central Locking Pedal	2pcs
✓ ABS Mattress Stopper	2pcs
✓ Hand Controller	1pc

Optional

<input type="checkbox"/> Manual CPR
<input type="checkbox"/> Electric CPR
<input type="checkbox"/> Bed Extension 20cm
<input type="checkbox"/> Back-up Battery
<input type="checkbox"/> Nurse Control Panel
<input type="checkbox"/> Emergency Stop Button
<input type="checkbox"/> Directional Wheels Function
<input type="checkbox"/> 5" and 6" Extra Side Rail For Full Protection
<input type="checkbox"/> Monkey Pole
<input type="checkbox"/> Balkan frame

Oasis Series

Y6w

Electric Bed
SKD-C

Technical Parameters

External Length	2130mm
External Width	1015mm
Bed Height Min.	365mm
Bed Height Max.	805mm
Mattress Platform	1925×880mm
Safe Working Load	250kg

Technical Configuration

✓ Linear Motor	3pcs
✓ 5" Covered Castors with Brake	4pcs
✓ Foldable Side Rail with ProErgo-Handle™	1 set
✓ ABS Oasis Bed Ends with Safe Lock	1 set
✓ PP Bed Platform, Removable & Washable	9pcs
✓ Ø70mm Wall Bumper Wheel	4pcs
✓ IV Pole 4 Hooks	1 set
✓ IV Pole Prevision	4pcs
✓ Drainage Hook	2pcs
✓ ABS Mattress Stopper	2pcs
✓ Hand Controller	1pc

Optional

<input type="checkbox"/> Manual CPR
<input type="checkbox"/> Electric CPR
<input type="checkbox"/> Bed Extension 20cm
<input type="checkbox"/> Back-up Battery
<input type="checkbox"/> Nurse Control Panel
<input type="checkbox"/> Emergency Stop Button
<input type="checkbox"/> Directional Wheels Function
<input type="checkbox"/> Monkey Pole
<input type="checkbox"/> Balkan frame



X9x

Electric ICU Bed SKD-C

Technical Parameters

External Length	2040-2240mm
External Width	1020mm
Bed Height Min.	500mm
Bed Height Max.	950mm
Mattress Platform	1820-2030×820mm
Safe Working Load	250kg

Optional

- ☐ Weighing System With LCD Screen
- ☐ Data Storage System
- ☐ Monkey Pole

Technical Configuration

✓ Unaki Motor	6pcs
✓ 6" Double Side Castors, Central Locking	4pcs
✓ ABS Side Rail with Integrated Control Panel	1set
✓ ABS Bed Ends One Step Removable	1set
✓ Ø75mm Wall Bumper Wheel	4pcs
✓ IV Pole 4 Hooks	1set
✓ IV Pole Prevision	4pcs
✓ Drainage Hook	4pcs
✓ 5" Directional Wheel Pedal	2pcs
✓ Central Locking Pedal	2pcs
✓ X-ray Translucent Mattress Platform	1set
✓ Back-up Battery	1set
✓ Night Light	1set
✓ Bed Push Handle for Emergency Transfer	1set
✓ Manual CPR	
✓ Electric CPR	



75°±10°

Back-rest



35°±10°

Knee-rest



290mm

Height Adjustment



12°

Reverse Trendelenburg



12°

Trendelenburg



CPR

Instant Release



0-90°

Angle Indicator



weight

Weighing System



Grace Care

D8d

Electric Bed
SKD-C

Technical Parameters

External Length	2280mm
External Width	1090mm
Bed Height Min.	435mm
Bed Height Max.	725mm
Mattress Platform	1925*900mm
Safe Working Load	250kg

Technical Configuration

✓ Linear Motor	4pcs
✓ 5" Double Side Castors, Central Locking	4pcs
✓ ABS Grace Side Rail with Integrated Control Panel	1set
✓ ABS Grace Bed Ends with Integrated LCD Control Panel	1set
✓ Ø85mm Wall Bumper Wheel	4pcs
✓ Trapeze Socket	2pcs
✓ IV Pole 4 Hooks	1set
✓ IV Pole Provision	4pcs
✓ Drainage Hook	2pcs
✓ Central Locking Pedal	2pcs
✓ Manual CPR	

Optional

<input type="checkbox"/> Electrical CPR
<input type="checkbox"/> Weighing System with LCD Screen
<input type="checkbox"/> Bed Extension
<input type="checkbox"/> Back-up Battery
<input type="checkbox"/> Cardiac Position
<input type="checkbox"/> ABS Bed Surface
<input type="checkbox"/> X-ray Bed Surface
<input type="checkbox"/> Monkey Pole
<input type="checkbox"/> Balkan Traction Frame
<input type="checkbox"/> Directional Wheel



SAIKANG



V8v

Electric Bed SKD-C



Technical Parameters

External Length	2085mm
External Width	1025mm
Bed Height Min.	455mm
Bed Height Max.	745mm
Mattress Platform	1925×900mm
Safe Working Load	250kg

Optional

- ☐ Battery Back-up
- ☐ Bed Extension
- ☐ Weighing System
- ☐ Hand Controller
- ☐ Night Light
- ☐ X-ray Translucent Mattress Platform (Back Section)
- ☐ Wheels Directional Function
- ☐ Utility Shelf
- ☐ Trapeze Bracket/Pole

Technical Configuration

✓ Linear Motor	4pcs
✓ 5" Double Side Castors, Central Locking	4pcs
✓ PP Side Rail Extra Height	1set
✓ PP Bed Ends with Safe Lock with Φ65 Bumper Wheels	1set
✓ PP Bed Platform with Integrated Mattress Retainer	1set
✓ Patent Angle Indicator For Back-rest And Trendelenburg Position	4pcs
✓ IV Pole 4 Hooks	1set
✓ IV Pole Prevision	4pcs
✓ Drainage Hook	2pcs
✓ Central Locking Pedal	2pcs
✓ Cardiac Chair Position	1pc
✓ Nurse Controller	1pc
✓ Integrated Control Panel	2set
✓ Manual CPR	
✓ Auto Break at 30°	
✓ Emergency Stop Button	

V8w

Electric Bed SKD-C



Technical Parameters

External Length	2085mm
External Width	1020mm
Bed Height Min.	455mm
Bed Height Max.	745mm
Mattress Platform	1925×900mm
Safe Working Load	250kg

Technical Configuration

✓ Linear Motor	4pcs
✓ Hand Controller	1pc
✓ 5" Handrail Castors	4pcs
✓ Foldable Side Rail with ProErgo-Handle™	1set
✓ PP Bed Ends with Safe Lock with Φ65 Bumper Wheels	1set
✓ PP Bed Platform with Integrated Mattress Retainer	1set
✓ IV Pole 4 Hooks	1set
✓ IV Pole Prewision	4pcs
✓ Drainage Hook	2pcs

Optional

<input type="checkbox"/> Battery Back-up
<input type="checkbox"/> Bed Extension
<input type="checkbox"/> Nurse Controller
<input type="checkbox"/> Emergency Stop Button
<input type="checkbox"/> Auto-regression Function (Back Section)
<input type="checkbox"/> Night Light
<input type="checkbox"/> Manual CPR
<input type="checkbox"/> X-ray Translucent Mattress Platform (Back Section)
<input type="checkbox"/> Utility Shelf
<input type="checkbox"/> Trapeze Bracket/Pole

V6v

Electric Bed SKD-C



Technical Parameters

External Length	2085mm
External Width	1025mm
Bed Height Min.	430mm
Bed Height Max.	740mm
Mattress Platform	1925×900mm
Safe Working Load	250kg

Technical Configuration


✓ Linear Motor	3pcs
✓ Hand Controller	1pc
✓ 5" Handrail Castors	4pcs
✓ PP Side Rail Extra Height	1set
✓ PP Bed Ends with Safe Lock with Φ65 Bumper Wheels	1set
✓ PP Bed Platform with Integrated Mattress Retainer	1set
✓ Patent Angle Indicator for Back-rest Position	2pcs
✓ IV Pole 4 Hooks	1set
✓ IV Pole Prewision	4pcs
✓ Drainage Hook	2pcs

Optional

<input type="checkbox"/> Battery Back-up
<input type="checkbox"/> Bed Extension
<input type="checkbox"/> Integrated Control Panel
<input type="checkbox"/> Nurse Controller
<input type="checkbox"/> Auto Break at 30°
<input type="checkbox"/> Emergency Stop Button
<input type="checkbox"/> Night Light
<input type="checkbox"/> Manual CPR
<input type="checkbox"/> X-ray Translucent Mattress Platform (Back Section)
<input type="checkbox"/> Utility Shelf
<input type="checkbox"/> Trapeze Bracket/Pole

V4w

Manual Bed
SK-CD1



Technical Parameters

External Length	2130mm
External Width	1020mm
Bed Height Min.	475mm
Bed Height Max.	725mm
Mattress Platform	1925×900mm
Safe Working Load	200kg

Technical Configuration


✓ Manual Cranks	4pcs
✓ 5" Hardshell Castors	4pcs
✓ Foldable Side Rail with ProErgo-Handle™	1 set
✓ PP Bed Ends with Safe Lock with Φ65 Bumper Wheels	1 set
✓ PP Bed Platform with Integrated Mattress Retainer	1 set
✓ IV Pole 4 Hooks	1 set
✓ IV Pole Provision	4pcs
✓ Drainage Hook	2pcs

Optional

<input type="checkbox"/> Auto-regression Function (Back Section)
<input type="checkbox"/> Utility Shelf
<input type="checkbox"/> Trapeze Bracket/Pole

V3v

Manual Bed
SK-CD1



Technical Parameters

External Length	2130mm
External Width	1025mm
Bed Height Min.	430mm
Bed Height Max.	740mm
Mattress Platform	1925×900mm
Safe Working Load	200kg

Technical Configuration

✓ Manual Cranks	3pcs
✓ 5" Hardshell Castors	4pcs
✓ PP Side Rail Extra Height	1 set
✓ PP Bed Ends with Safe Lock with Φ65 Bumper Wheels	1 set
✓ PP Bed Platform with Integrated Mattress Retainer	1 set
✓ Patent Angle Indicator For Back-rest Position	2pcs
✓ IV Pole 4 Hooks	1 set
✓ IV Pole Provision	4pcs
✓ Drainage Hook	2pcs

Optional

<input type="checkbox"/> Utility Shelf
<input type="checkbox"/> Trapeze Bracket/Pole



V6k

Electric Bed SKD-C

RAL 110 715 96 RAL 5003 RAL 6005 RAL 3016 Pear Oak White Oak Red Oak

Technical Parameters

External Length	2085mm
External Width	995mm
Bed Height Min.	420mm
Bed Height Max.	730mm
Mattress Platform	1925×900mm
Safe Working Load	200kg

Technical Configuration

✓ Linear Motor	3pcs
✓ Hand Controller	1pc
✓ 5" Covered Castors	4pcs
✓ Foldable Side Rail	1set
✓ PP Bed Ends with Safe Lock with Φ65 Bumper Wheels	1set
✓ PP Bed Platform with Integrated Mattress Retainer	1set
✓ IV Pole 4 Hooks	1set
✓ IV Pole Prevision	4pcs
✓ Drainage Hook	2pcs

Optional

<input type="checkbox"/> Battery Back-up
<input type="checkbox"/> Bed Extension
<input type="checkbox"/> Nurse Controller
<input type="checkbox"/> Emergency Stop Button
<input type="checkbox"/> Auto-regression Function (Back Section)
<input type="checkbox"/> Night Light
<input type="checkbox"/> Manual CPR
<input type="checkbox"/> X-ray Translucent Mattress Platform (Back Section)
<input type="checkbox"/> Utility Shelf
<input type="checkbox"/> Trapeze Bracket/Pole



360° View

V3k

Manual Bed
SK-CD1

Technical Parameters

External Length	2130mm
External Width	995mm
Bed Height Min.	420mm
Bed Height Max.	730mm
Mattress Platform	1925×900mm
Safe Working Load	200kg

Technical Configuration

✓ Manual Crank	3pcs
✓ 5" Covered Castors	4pcs
✓ Foldable Side Rail	1 set
✓ PP Bed Ends with Safe Lock with Φ65 Bumper Wheels	1 set
✓ PP Bed Platform with Integrated Mattress Retainer	1 set
✓ IV Pole 4 Hooks	1 set
✓ IV Pole Prevision	4pcs
✓ Drainage Hook	2pcs

Optional

<input type="checkbox"/> Auto-regression Function (Back Section)
<input type="checkbox"/> Utility Shelf
<input type="checkbox"/> Trapeze Bracket/Pole

V2k

Manual Bed
SK-C1

Technical Parameters

External Length	2130mm
External Width	995mm
Bed Height	500mm
Mattress Platform	1925×900mm
Safe Working Load	200kg

Technical Configuration

✓ Manual Crank	2pcs
✓ 5" Covered Castors	4pcs
✓ Foldable Side Rail	1 set
✓ PP Bed Ends with Safe Lock with Φ65 Bumper Wheels	1 set
✓ PP Bed Platform with Integrated Mattress Retainer	1 set
✓ IV Pole 4 Hooks	1 set
✓ IV Pole Prevision	4pcs
✓ Drainage Hook	2pcs

Optional

<input type="checkbox"/> Auto-regression Function (Back Section)
<input type="checkbox"/> Utility Shelf
<input type="checkbox"/> Trapeze Bracket/Pole



User's Safety and Comfort from Long Term Care Bed

The physical health and psychological well-being are hot topics for elderly people. Safety is an important quality feature for a nursing bed. Beds' systematic fall prevention can reduce the risk of falling out of the bed.

Our innovative solutions can reduce falls as well as the consequences associated with falls. It is a wonderful match of safety, mobility, height adjustable and coziness. All these features cater to the needs of people who are spend a lot time in bed. The adjustability makes beds more appropriate than a standard bed for providing comfortable support and enabling caregivers to deliver quality care.

Beech Wood Cherry Wood Oak Oak Dark Oak

HB421

Electric Homecare bed
SKD-C
Size: 2150×920×185-762mm



HB422

Electric Homecare bed
SKD-C
Size: 2150×1000×200-745mm



HB425

Electric Homecare bed
SKD-C
Size: 2150×1000×185-762mm




HB426

Electric Homecare bed
SKD-C
Size: 2200×1050×345-927mm





100% View



HC300 is specially made for home care and long term care environments and offers superior comfort and usability features for its users.



HC300

Electric Homecare bed
SKD-C

Size: 2037×996×250-580mm



HC200

Electric Homecare bed
SKD-C

Size: 2037×996×250-580mm



10-piece steel plates

The HC nursing beds are designed with light weight, practical structure. The bed deck is composed of 10-piece high strength steel plates. And with ventilator holes which can improve air flow.



Head & Foot Board

The headboard & footboard are made of MDF, which meets the environmental protection requirements of E1 level. HF board edges are round, designed for easy grasping, additional leather cover is optional with different colors to complement customer's decoration and furnishings.



Smart slow start and stop

Smart slow start and stop: Running speed of the back plate is reduced to half, at the beginning and ending of the movement, so as to reduce the sudden impact on the elderly.



One-button function

The specific one-button function makes the bed more convenient to operate. The remote control is equipped with a relaxation position, reading position, out of bed position.



GB4w

Manual Orthopedic Traction Bed

SKG-4C



Technical Parameters

External Length	2180mm
External Width	960mm
Bed Height Min.	460mm
Bed Height Max.	710mm
Matress Platform	1925×830mm

Technical Configuration

✓ Linear Motor	3pcs	✓ IV Pole	1 pc
✓ 5" Hardshell Castors	4pcs	✓ IV Pole Prevision	4pcs
✓ Aluminium Alloy Foldable Side Rail	1 set	✓ Drainage Hook	2pcs
✓ Wood Bed Ends	1 set	✓ ABS Mattress Holder	2pcs
✓ Wood Bed Side Panel	1 set	✓ Aluminium Balkan Frame	1 set



Design

Use new and strong inserted ABS fire-proof panel avoid off. Easy to be cleaned and installed.



Steels | Painting

Baosteel from Fortune Global 500. 11 Process epoxy painting, ASTM testing anti-bacterial, paint thickness 0.12mm, brightness 60°, paint can resist 50kg impact.



Welding | Chassis

Panasonic robotic ensure, 360° full smooth wedding Solid bed base made of 5.0 mm stamped rotating plate.



Hand Crank

Retractable hand crank, 100,000 times remaining intact. Extreme protection device, sounds when tipped to the max position.



Balkan Frame

Aluminum Balkan Frame, which can adjust the longitudinal, lateral and height positions in all directions, and can accurately adjust the most comfortable posture of the patient.



Head & Foot Board

Full ABS material, head and foot board can be reversible.

CX2x

Children Bed

SK-C2



Technical Parameters

External Length	1880mm
External Width	900mm
Bed Height	585mm
External Height	1260mm
Mattress Platform	1770×720mm

Technical Configuration

✓ Manual Cranks	2pcs	✓ IV Pole	1pc
✓ 5" Covered Castors	4pcs	✓ Drainage Hook	2pcs
✓ Full Length Side Rail	1set	✓ Buffer Air Spring	4pcs
✓ Plastic Bed Ends	1set		
✓ IV Pole Prevision	4pcs		

CR0q

Children Bed

SK-A2



Technical Parameters

External Length	1600mm
External Width	740mm
Bed Height	500mm
External Height	1040mm
Mattress Platform	1580×660mm

Technical Configuration

✓ 5" Covered Castors	4pcs	✓ Drainage Hook	2pcs
✓ Full Length Side Rail	1set		
✓ Integrated Bed Ends	1set		
✓ IV Pole	1pc		
✓ IV Pole Prevision	4pcs		

CR2q

Children Bed

SK-C2



Technical Parameters

External Length	1820mm
External Width	860mm
Bed Height	600mm
External Height	1100mm
Mattress Platform	1600×680mm

Technical Configuration

✓ Manual Cranks	2pcs	✓ IV Pole Prevision	4pcs
✓ Full Length Side Rail	1set	✓ Drainage Hook	2pcs
✓ Integrated Bed Ends	1set	✓ Shoes Holder	1pc
✓ Anti-slip Bed Feet Cover	4pcs		
✓ IV Pole	1pc		

CD3q

Children Bed


SK-CD2



CT8k

Electric Children Bed

SKD-C



Technical Parameters

External Length	1835mm
External Width	780mm
Bed Height Min.	480mm
Bed Height Max.	770mm
Mattress Platform	1625×650mm

Technical Configuration

✓ Linear Motor	4pcs	✓ IV Pole	1pc
✓ 5" Double Side Castors, Central Locking	4pcs	✓ Drainage Hook	2pcs
✓ Aluminium Alloy Side Rail with Extra Protection	1set		
✓ Aluminium Alloy Bed Ends	1set		
✓ IV Pole Prevision	4pcs		

X05-3

Children Bed

SK-C2



Technical Parameters

External Length	1800mm
External Width	800mm
Bed Height	700mm
External Height	1100mm
Mattress Platform	1435×690mm

Technical Configuration

✓ Manual Cranks	2pcs	✓ IV Pole	1pc
✓ 5" Covered Castors	4pcs	✓ Drainage Hook	2pcs
✓ Stainless Steel Full Length Side Rail	1set		
✓ Stainless Steel Bed Ends	1set		
✓ IV Pole Prevision	4pcs		

Technical Parameters

External Length	2030mm
External Width	1045mm
Bed Height Min.	500mm
Bed Height Max.	770mm
Mattress Platform	1705×900mm

Technical Configuration

✓ Manual Cranks	3pcs	✓ IV Pole	1pc
✓ 5" Hardshell Castors	4pcs	✓ Drainage Hook	2pcs
✓ Aluminium Alloy Full Length Side Rail	1set		
✓ ABS Groce Bed Ends with Safe Lock	1set		
✓ IV Pole Prevision	4pcs		



X01-1

Baby Crib
SK-AD3



Technical Parameters

External Length	875mm
External Width	525mm
Bed Height Min.	650mm
Bed Height Max.	900mm

Technical Configuration

✓ Hydraulic Pump	1pc
✓ Bassinet	1pc
✓ 3" Double Side Castors	4pcs

X01-5

Baby Crib
SK-A3



X01
Stainless Steel



Powder Coated Steel

Technical Parameters

External Length	810mm
External Width	520mm
Bed Height	990mm

Technical Configuration

✓ Bassinet	1pc	✓ IV Pole	1set
✓ 3" Covered Castors	4pcs	✓ Utility Basket	1pc
✓ IV Pole Prevision	2pcs	✓ Adjustable Air Spring	1pc

X20

Pediatric Exam Table
SKZ-A1

Technical Parameters

External Length	1800mm
External Width	650mm
Bed Height	830mm



Z0n

Electric Bath Bed
SKD-A



Z0y
Hydraulic



Technical Parameters

External Length	1950mm
External Width	780mm
Bed Height Min.	600mm
Bed Height Max.	830mm


Technical Configuration

✓ Linear Motor	1pc
✓ 3" Covered Castors	4pcs
✓ Stainless Steel Side Rail	1set
✓ Stainless Steel Bed Ends	1set

X12

Electric Examination Table
SKD-A

230mm
Height
Adjustment
500-730mm



Technical Parameters

External Length	1900mm
External Width	600mm
Table Height Min.	500mm
Table Height Max.	730mm

Technical Configuration

✓ Linear Motor	1pc
✓ Anti-slip Bed Feet Cover	4pcs
✓ Removable Pillow	1pc

X14

Hydraulic Examination Table
SKZ-CD2

75°
Back-rest



400mm
Height
Adjustment
500-900mm



Technical Parameters

External Length	1940mm
External Width	700mm
Table Height Min.	490mm
Table Height Max.	890mm

Technical Configuration

✓ Hydraulic Pump	1pc
✓ Adjustable Air Spring	2pcs
✓ 5" Covered Castors	4pcs
✓ Hydraulic pedal	2sets

X26

Electric Examination Table
SKD-A

75°
Back-rest



370mm
Height
Adjustment
480-850mm



Technical Parameters

External Length	1810mm
External Width	680mm
Table Height Min.	480mm
Table Height Max.	850mm

Technical Configuration

✓ Pillow	1pc	✓ Hand Controller	1pc
✓ Motor	1pc	✓ Paper Roll Holder	1pc
✓ Mechanical Lever to Adjust Back Rest	1pc		
✓ 2" Casters	2pcs		

X36

Electric Examination Table
SKD-A

75°/54°
Back Plate
up/down



33°/32°
Leg Plate
up/down



450mm
Height
Adjustment
460-910mm



Technical Parameters

External Length	1890mm
External Width	860mm
Table Height Min.	460mm
Table Height Max.	910mm

Technical Configuration

✓ Pillow	1pc	✓ Brake Pedal	4pcs
✓ Motor	1pc	✓ Hand Controller	1pc
✓ Gas spring	2pcs	✓ Manual Operator	1pc
✓ 3" Casters	4pcs		

X07a
Examination Table
SKZ-A1



Technical Parameters

External Length	1900mm
External Width	600mm
Table Height	680mm

Technical Configuration

✓ Anti-slip Bed Feet Cover	4pcs

X09
Examination Table
SKZ-B1




Technical Parameters

External Length	1900mm
External Width	600mm
Table Height	680mm

Technical Configuration

✓ Anti-slip Bed Feet Cover	4pcs
✓ Mechanical Lever to Adjust Back Rest	1set

A048
Gynecological Examination Table
SKF-B



Technical Parameters

External Length	1810mm
External Width	690mm
Table Height	880mm

Technical Configuration

✓ Drawer	4pcs	✓ Lithotomy Rod	1set
✓ Step Stool	1pc	✓ Stainless Steel Basin	1pc
✓ Retractable Foot Board	1pc	✓ Paper Roll Holder	1pc

A046

Gynecological Table
SKF-B


Technical Parameters

External Length	1900mm
External Width	600mm
Table Height	750mm

Technical Configuration

✓ Crank Wheel	2pcs	✓ Stirrup	1set	✓ Anti-slip Bed Feet Cover	4pcs
✓ Lithotomy Rod	1set	✓ Hand Bar	1set	✓ Mechanical Lever	2pcs
✓ Shoulder Holder	1set	✓ Arm Holder	1set	✓ 3" Castors with Brake	4pcs

A045-1

Gynecological Table
SKF-A

Technical Parameters

External Length	1820mm
External Width	590mm
Table Height	800mm

Technical Configuration

✓ Lithotomy Rod	1set	✓ Mechanical Lever	2sets
✓ Hand Bar	1set	✓ Stainless Steel Basin	1pc
✓ Anti-slip Bed Feet Cover	4pcs		

A99-8

Electric Gynecological Exam Couch
SKL-B

Technical Parameters

External Length	1400mm
External Width	740mm
Table Height Min.	610mm
Table Height Max.	1100mm

Technical Configuration

✓ Hand Bar	2pcs
✓ Leg Support with Cushion	2pcs
✓ Foot Pedal	1pc
✓ Stainless Steel Basin	1pc

A99-6

Electric Gynecological Exam Couch
SKL-B

Technical Parameters

External Length	1200mm
External Width	500mm
Table Height Min.	500mm
Table Height Max.	1000mm

Technical Configuration

✓ Simple Side Rail	1set
✓ Leg Holder	1set
✓ 5" Castors 4	4pcs
✓ Foot Pedal	1pc



A99-7
Electric Obstetric Bed
SKL-A

250mm
Height
Adjustment

55°/18°
Back Plate
up/down

35°/5°
Seat Plate
up/down

90°
Leg Plate
Swing Out

ELECTRIC OBSTETRIC BED

Frame

Base cover made of high quality ABS. It can be sterilized at the temperature above 100 °C. The tensile strength up to 30 MPa.

Base structure is made of high-quality seamless steel, robot welded with high strength.

Battery

Battery could support about 80-120 complete operation cycles. Working time could be last for one week according to useage and help to keep the operating table running in the event of power failure.



Bed Width

Seamless foam mattress, with waterproof PU cover. Easy cleaning and disinfection.

Arm Plate

Flexible arm plate adjustment.

Hand Controller

Electrical braking system. Controlled by hand controller.

Technical Parameters

External Length	1980mm
External Width	720mm
Table Height Min.	630mm
Table Height Max.	880mm

Technical Configuration

✓ Bed Frame	1pc	✓ Hand Bars	2pcs	✓ Caster	1pc
✓ Foot Pedal	1pc	✓ Leg Support with Cushion	2pcs		
✓ Hand Controller	1pc	✓ Arm Support with Cushion	2pcs		
✓ Anesthesia Screen	1pc	✓ Stainless Steel Basin	1pc		

A99-5

Electric Obstetric Bed
SKL-A

Technical Parameters

External Length	1880mm
External Width	660mm
Table Height Min.	680mm
Table Height Max.	940mm

Technical Configuration

✓ Arm Support	2pcs	✓ Tightener	2pcs	✓ Slider	4pcs
✓ Handle	2pcs	✓ Hand Controller	1pc		
✓ Leg Support	2pcs	✓ Foot Switch	1pc		
✓ Stirrup	2pcs	✓ Stainless Steel Basin	1pc		

A99-13

Electric Obstetric Bed
SKL-A

Technical Parameters

External Length	2100mm
External Width	600mm
Table Height Min.	750mm
Table Height Max.	1000mm

Technical Configuration

✓ Armrest	2pcs	✓ SST Rinsing Basin	1pc
✓ Knee Cushion	2pcs	✓ Pedal Controller	1pc
✓ Pull Handle	2pcs	✓ LED Examination Light	1pc
✓ Anesthesia Frame	1pc		

A105

Manual Obstetric Table
SKF-C

Technical Parameters

External Length	1900mm
External Width	600mm
Table Height Min.	700mm
Table Height Max.	900mm

Technical Configuration

✓ Bed Frame	1pc	✓ Leg Support with Cushion	2pcs
✓ Mattress	1pc	✓ Arm Support with Cushion	2pcs
✓ Anesthesia Screen Frame	1pc	✓ Stainless Steel Basin	1pc
✓ Hand Bars	2pcs		

A98-4

Electric Obstetric Bed
SKL-A

Technical Parameters

External Length	1920mm
External Width	920mm
Table Height Min.	650mm
Table Height Max.	900mm

Technical Configuration

✓ Mattress	1pc	✓ Slider	2pcs	✓ Castors	4pcs
✓ Handle	2pcs	✓ Tightener	2pcs	✓ Central Locking Pedal	2pcs
✓ Leg Support	2pcs	✓ Hand Controller	1pc		
✓ ABS Side Rail	2pcs	✓ Stainless Steel Basin	1pc		



A501
Electric-Hydraulic
Operating Table
SKT-B

A501 is an advanced, multi-functional electric hydraulic operating table. Its table top is made of high-strength synthetic board and can transmit X-rays. In addition, the platform can be equipped with a 180° rotation function to realize head and leg interchange, which is applicable for C-arm perspective and placement of some special positions.

15° Reverse Trendelenburg	0° Trendelenburg	15° Lateral Tilt left/right	40°/90° Head Plate up/down	85°/40° Back Plate up/down	15°/90° Leg Plate up/down
80mm Kidney Bridge Up	360mm Horizontal Sliding Crank Arm	180° Horizontal Rotation	95° Flex	220° Reflex	



Technical Parameters

External Length	2125mm
External Width	550mm
Table Height Min.	600mm
Table Height Max.	950mm

Technical Configuration

✓ Bed Frame	1pc	✓ Arm Holder with Cushion	2pcs	✓ Kidney Bridge Handle	1pc
✓ Mattress	1set	✓ Leg Support with Cushion	2pcs		
✓ Anesthesia Screen Frame	1pc	✓ Foot Support with Cushion	2pcs		
✓ Shoulder Holder with Cushion	2pcs	✓ Hand Controller	1pc		



A106-2
Electric-Hydraulic Operating Table
SKT-B

A106-2 is a highly automated electric-hydraulic operating table, which can meet the needs of various surgical operations. The unique design of a super-low position of 500mm, aimed for the brain surgery and ophthalmic surgery.

25°
Reverse
Trendelenburg

25°
Trendelenburg

20°
Lateral Tilt
left/right

45°/90°
Head Plate
up/down

80°/20°
Back Plate
up/down

15°/90°
Leg Plate
up/down

90°
Foot Plate
Horizontal
Positioning

100mm
Kidney
Bridge Up



Technical Parameters

External Length	2080mm
External Width	550mm
Table Height Min.	500mm
Table Height Max.	800mm

Technical Configuration

√ Anaesthesia Screen	1pc	√ Foot Support	1pc	√ Fixing Clamp	7pcs
√ Hand Remote	1pc	√ Mattress	1set	√ Power Line	1pc
√ Kidney Bridge Handle	1pc	√ Shoulder Support	1pair		
√ Arm Support	1pc	√ Long Fixing Clamp	1pair		



360° View

A302

Electric Operating Table SKL-A

A302 is a highly automated electric Multi-functional operating table. Its elevation, lowering, Trendelenburg, reverse Trendelenburg, lateral tilt, back folding, vertical translation and leg folding are all electrically operated. The operating table is equipped with a digital display hand controller to make operation easy and intuitive. With imported large diameter castors, smooth movement, stable brake system. The standard dual control system guarantees the safety of operation.



Material

The base material for the table is high thickness 304 stainless steel, seamless welding by robot, with modern appearance and more stable.



Mattress

Polyurethane foam mattress, waterproof and seamless, avoid liquid in. The connection between table and mattress is buckle joint, easy to disassemble and install, easy to clean.



Head Plate

Foldable head plate, support for more type of surgical operations.



Leg Plate

The up and down folding actions of the leg plate are all completed by electric. The leg plate is removable, effortless and smooth.



Hand Controller

All functions controlled by both hand controller or control panel on table base column. With position LCD display on hand controller, easy to operate.



Battery Could Support

Battery could support about 80-120 complete operation cycles. Working time could be last for one week according to usage and help to keep the operating table tuning in the event of power failure.

Technical Parameters

External Length	2150mm
External Width	555mm
Table Height Min.	700mm
Table Height Max.	1000mm

Technical Configuration

✓ Bed Frame	1 pc	✓ Arm Holder	2 pcs	✓ Sliding Block	8 pcs
✓ Mattress	1 set	✓ Leg Support	2 pcs	✓ Spanner	5 pcs
✓ Anesthesia Screen Frame	1 pc	✓ Foot Support	2 pcs		
✓ Shoulder Holder	2 pcs	✓ Hand Controller	1 pc		

A301

Electric Operating Table
SKL-A

A301 is a multi-functional electric operating table. Its mattress is made of polyurethane molded material, can be sterilized and cleaned at high temperature; and there are no gaps so to avoid bacterial growth. The platform is made of high-purity phenolic resin board applicable for C-arm, and can be sided up to 350mm in one direction. The X-ray photography can carry out the whole body radiography examination without moving the patient, no blind spots or shadows.



Technical Parameters

External Length	2150mm
External Width	555mm
Table Height Min.	750mm
Table Height Max.	1000mm

Technical Configuration

✓ Bed Frame	1pc	✓ Arm Holder	2pcs	✓ Sliding Block	8pcs
✓ Mattress	1set	✓ Leg Support	2pcs	✓ Spinner	5pcs
✓ Anesthesia Screen Frame	1pc	✓ Foot Support	2pcs		
✓ Shoulder Holder	2pcs	✓ Hand Controller	1pc		

A100-4

Electric Operation Table
SKL-A

A100-4 is an electric operating table that can meet the needs of various surgeries. It is driven by low-voltage DC motor, with safe and stable performance. Super wide table top, super long sliding distance. It can be applicable for C-arm. It has a high level of automation, low noise, strong reliability, and imported original parts for key parts. It is a perfect and ideal operating table.



Technical Parameters

External Length	2070mm
External Width	550mm
Table Height Min.	700mm
Table Height Max.	1000mm

Technical Configuration

✓ Bed Frame	1pc	✓ Arm Holder with Cushion	2pcs
✓ Mattress	1set	✓ Leg Support with Cushion	2pcs
✓ Anesthesia Screen Frame	1pc	✓ Foot Support with Cushion	2pcs
✓ Shoulder Holder with Cushion	2pcs	✓ Kidney Bridge Handle	1pc

A100-4A

Electric Operating Table

SKL-A

A100-4A is an electric operating table with super load bearing and extra-long sliding distance, which can be applicable for C-arm. The built-in large-capacity battery can continue to support all operations in the event of a power failure. And the operating table adopts a dual control system to ensure the safety of the operation.



Technical Parameters

External Length	2100mm
External Width	550mm
Table Height Min.	700mm
Table Height Max.	1000mm

Technical Configuration

Bed Frame	1pc	Arm Holder with Cushion	2pcs	Kidney Bridge Handle	1pc
Mattress	1set	Leg Support with Cushion	2pcs		
Anesthesia Screen Frame	1pc	Foot Support with Cushion	2pcs		
Shoulder Holder with Cushion	2pcs	Hand Controller	1pc		

A2000

Electric Operating Table

SKL-A

A2000 is a basic electric operating table with a very high cost effective. It can fulfill the requirement of general surgeries. Its built-in battery can continue to support all operations in the event of a power failure.



Technical Parameters

External Length	2100mm
External Width	480mm
Table Height Min.	680mm
Table Height Max.	930mm

Technical Configuration

Bed Frame	1pc	Arm Holder	2pcs	Sliding Block	10pcs
Mattress	1set	Leg Support	2pcs	Spanner	4pcs
Anesthesia Screen Frame	1pc	Foot Support	2pcs		
Shoulder Holder	2pcs	Hand Controller	1pc		

A205

Manual Operating Table
SKF-C

Technical Parameters

External Length	2000mm
External Width	480mm
Table Height Min	720mm
Table Height Max	960mm

Technical Configuration

✓ Arm support	2 pcs	✓ Anesthetic Screen	1 pc
✓ Leg support	2 pcs		
✓ Body Support	2 pcs		
✓ Shoulder Support	2 pcs		

A3001

Manual Operating Table
SKF-C

Technical Parameters

External Length	2000mm
External Width	480mm
Table Height Min	790mm
Table Height Max	1040mm

Technical Configuration

✓ Anesthesia Screen Frame	1 pc	✓ Leg Support with Cushion	2 pcs
✓ Mattress	1 set	✓ Shoulder Support with Cushion	2 pcs
✓ Kidney Bridge Handle	1 pc	✓ Body Support with Cushion	2 pcs
✓ Arm Support with Cushion	2 pcs		

Multi-functional Orthopedic Frame



Gel Mattress

Connecting Multi-functional
Orthopedic FrameCarbon Fiber Multi-functional
Orthopedic Frame

Multi-purpose Head Rack



Shoulder Traction Frame



Stirrup Pneumatic Power Leg Rack



Bow Spine Operating Supporter



Mechanical Power IV Pole



Operating Table Accessories Cart



Knee Surgery Rack



Double Table Top



Prone Body Supporting Frame



Arm Microsurgical Rack



Bed Surface Widening Module



External Kidney Bridge



Dirt Basin



Body Safe Strap



Mattress [OEM Size]

SKP011



SKP104



SKP008



SKP006

Air Mattress



Mattress

Arrangement	2.5" bubble pad	Dimension	2000*900*70mm
Material	medical grade PVC	Weight	3kg

SKP013

Air Mattress



Mattress

Standard Size	80"×32"×5"	Bottom	tarpaulin /nylon PU
Top cover	Nylon PU /strength	Mattress weight	5.5-7kg
Air cells	Nylon / PVC/TPU	Weight Capacity	135-145kg

SKP012

Air Mattress



Mattress

Standard Size	80"×32"×6"	Bottom	tarpaulin /nylon PU
Top cover	Nylon PU /strength	Mattress weight	5.5-6.5kg
Air cells	Nylon / PVC/TPU	Weight Capacity	135-145kg

SKP015

ICU Air Mattress



Mattress

Standard Size	80"×32"×8"	Bottom	tarpaulin /nylon PU
Top cover	Nylon PU /strength	Mattress weight	10kg
Air cells	Nylon / PVC/TPU	Weight Capacity	180kg

SKS002

ABS Bedside Table

480×465×745mm(L×W×H)

**SKS002-W**

ABS Bedside Table

480×465×785mm(L×W×H)

**SKS002-5**

ABS Bedside Table

480×480×760mm(L×W×H)

**SKS003**

ABS Bedside Table

480×460×900mm(L×W×H)

**SKS005**

Steel Bedside Table

480×400×790mm(L×W×H)

**SKS007-3**

Steel Bedside Table

440×460×810mm(L×W×H)

**SKS007-4**

Steel Bedside Table

440×460×810mm(L×W×H)

**SKS036**

Bedside Table

440×440×790mm(L×W×H)

**SKS008**

ABS+Steel Bedside Table

500×510×705mm(L×W×H)

**SKS013-1**

ABS+Steel Bedside Table

575×500×750mm(L×W×H)

**SKS006-1**

Stainless Steel Bedside Table

475×390×760mm(L×W×H)

**SKS020**

Solid Wood Bedside Table

500×430×650mm(L×W×H)

**SKS025**

Phenolic Resin Bedside Table

475×475×800mm(L×W×H)

SKH046-11

Overbed Table

Technical Parameters

800×395×870-1170mm

Technical Configuration

✓ Air Spring for Height Adjustment

✓ Phenolic Resin Table Top

300mm
Height
Adjustment**SKH046-2**

Overbed Table

400mm
Height
Adjustment

Technical Parameters

890×440×765-1050mm

Technical Configuration

✓ Air Spring for Height Adjustment

✓ ABS Table Top with Cup Socket

**SKH046-14**

Overbed Table

300mm
Height
Adjustment

Technical Parameters

890×400×890-1190mm

Technical Configuration

✓ Air Spring for Height Adjustment

✓ ABS Table Top with Cup Socket and Push Handle

**SKH042**

Overbed Table

400mm
Height
Adjustment

Technical Parameters

800×410×730-1130mm

Technical Configuration

✓ Air Spring for Height Adjustment

✓ MDF Table Top

**SKH046-1**

Telescopic Dining Board

150mm
Telescopic

Technical Parameters

870-1020×350×28mm

Technical Configuration

✓ Hook



SKE001

Attendant Chair

Technical Parameters

Seat Size	680×720×940mm
Bed Size	680×1900×590mm

Technical Configuration

✓ Arm	1set
✓ 2" Castors	4pcs
✓ 3" Castors	2pcs
✓ Anti-slip Feet Cover	4pcs
✓ Pillow	1pc

**SKE001-3**

Luxury Attendant Chair

Technical Parameters

Seat Size	710×720×950mm
Bed Size	710×1940×420mm

Technical Configuration

✓ Wood Arm Frame	1set
✓ 2" Castors	2pcs
✓ Anti-slip Feet Cover	2pcs
✓ Pillow	1pc

**SKE001-4**

Luxury Attendant Sofa

Technical Parameters

Seat Size	1350×900×650mm
Bed Size	1900×1350×500mm

**SKE014-5**

Luxury Attendant Sofa

Technical Parameters

Seat Size	900×900×650mm
Bed Size	1900×900×500mm

**SKE-135**

Electric Two Function Dialysis Chair

0°-65°
Back Section
Adjustment10°-90°
Leg Section
Adjustment

Technical Parameters

1830×580×560mm

Technical Configuration

- ✓ Hand Controller
- ✓ Power Line
- ✓ 3" Silent Caster
- Small Dining Table

**SKE-131**

Electric Two Function Dialysis Chair

0°-70°
Back Section
Adjustment0°-60°
Leg Section
Adjustment

Technical Parameters

1880×580×560mm

Technical Configuration

- ✓ Hand Controller
- ✓ Power Line
- ✓ 3" Silent Caster
- Small Dining Table

**SKE-133**

Electric Three Function Dialysis Chair

-8°-70°
Back Section
Adjustment0°-75°
Leg Section
Adjustment550-750mm
Height
Adjustable

Technical Parameters

1800×580×(550-750)mm

Technical Configuration

- ✓ Hand Controller
- ✓ Power Line
- ✓ 3" Silent Caster
- ✓ With CPR
- Small Dining Table

**SKE-134**

Electric Three Function Dialysis Chair

0°-70°
Back Section
Adjustment20°-75°
Leg Section
Adjustment540-640mm
Height
Adjustable

Technical Parameters

1950×580×(540-640)mm

Technical Configuration

- ✓ Hand Controller
- ✓ Power Line
- ✓ 3" Silent Caster
- Small Dining Table



SKE-132

Electric Dialysis Chair



Technical Parameters

1900×700×550mm

Technical Configuration

- ✓ Hand Controller
- ✓ 3" Silent Caster
- Small Diving Table



SKE087

Reclining Chair

Technical Parameters

Seat Size: 820×850×1070mm
Bed Size: 820×1600×450mm



Teabag Mat



SKE091

Blood Donation Chair

Functions

Back-rest: 0°-70°
Leg Rest Down Fold: 0°-90°

Technical Parameters

Seat Size: 560×1030×1200mm
Bed Size: 560×1650×470mm

Technical Configuration

- ✓ Arm: 2pcs
- ✓ Anti-slip Feet Cover: 4pcs



SKE090

Blood Donation Chair

Technical Parameters

Seat Size: 910×480×960mm

Technical Configuration

- ✓ Arm: 2pcs
- ✓ Anti-slip Feet Cover: 4pcs



SKE004-1

Transfusion Chair

Technical Parameters

Seat Size: 670×800×1050mm

Technical Configuration

- ✓ Anti-slip Feet Cover: 4pcs
- ✓ Wooden Armrest: 2pcs
- ✓ IV Pole Previsions: 2pcs
- ✓ IV Pole: 1 set



SKE005

Transfusion Chair

Technical Parameters

Seat Size: 750×1000×1050mm
Bed Size: 750×1550×450mm

Technical Configuration

- ✓ Anti-slip Feet Cover: 4pcs
- ✓ Wooden Armrest: 2pcs
- ✓ IV Pole Previsions: 2pcs
- ✓ IV Pole: 1 set
- ✓ Utility Basket: 1 pc



Treat-waiting Chair



SKE008-1

Material: Powder Coated
1780×740×800mm(L×W×H)



SKE006-1

Material: PVC
1655×675×870mm(L×W×H)



SKE008

Material: Chromed Steel
1800×740×800mm(L×W×H)



SKE006-2

Material: PU
1765×675×905mm(L×W×H)



SKE901

Pediatric Infusion Chair



INFUSION ROOM

Preparation of infusion chair or infusion bed depends on the patient's condition. The selection of material has the characteristics of environmental protection, durability, alkali resistance and disinfectant cleaning. It is a smooth metal surface that scratches the patient.

SKE902

Infusion Sofa



SKE903

Infusion Sofa



SKE904

Pediatric Infusion Chair



SKE905

Pediatric Nebulization Table





360° View

SKZ412
Office Table
Size: 1800×800×750mm



SKZ412-1
Side Cabinet
Size: 800×400×650mm



SKZ413
Meeting Table
Size: Ø800×750mm



SKH091
Cabinet
Size: 900×400×1800mm



SKE702
Office Chair
Mesh Cloth with Nylon Base
Size: 480×430×460-540mm



SKE055-8
Office Chair
Mesh Chair Back high-quality Mesh with
Recycled Sponge, Plastic Handrail
Size: 480×470×900mm



SKE020-1
Step Stool
Size: 490×405×170/340mm



SKH049
Stainless Steel Screen
Size: 2000×1870mm



SK-CL001 **SK-CL501**
Curtain Track



SKH041-3
I.V. Stand
Size: 500×490×1350-2250mm



SKH041(1)
I.V. Stand
Size: Ø520×1350-2250mm



SKH034
Trash Can



12L

SKH034-5
Kick Bucket
Size: Ø400×H360×14L



14L

SKE180
Electric Wheelchair

Size: 1150×450×950mm

**SKE030**
WheelchairPowder Coated Steel Structure
Seat Dimension: 295×455mm
Size: 1035×450×670mm**SKE031**
Commode ChairHeight Adjustable, with Foot Rest
Aluminium Folding Frame, Plastic
Size: 880×510×300mm**SKE053-2**
Training ChairWriting Board (370×245mm)
Size: 560×450×750mm**SKE052**
Office ChairPowder Coated Steel Structure
Back and Seat Spun with the Elastic Fabric
Size: 495×470×720mm**SKE711**
Office ChairMesh Chair Back, high-quality Mesh with Recycled Springs, Plastic Handrail
Size: 520×450×920, 680mm**SKE712**
Office ChairMesh Chair Back, high-quality Mesh with Recycled Springs, plastic Handrail
Size: 550×450×1000mm**SKE802**
Dental ChairChrome Plated Structure with Polyurethane Foam Seat
Size: 6750×430×160mm**SKE015-2**
Nurse StoolChrome Plated Structure with PU Foam Seat Top
Size: 4340×490×680mm**SKE801**
Nurse StoolStainless Steel Base with PU Leather
High Density Sponge Seat
Size: 6710×490×630mm**SKE015-5**
Nurse StoolChrome Plated Structure with PU Leather
High Density Sponge Seat
Size: 6710×445×575mm**SKE013**
Nurse StoolChrome Plated Structure with PU Foam Seat Top
Size: 6710×445×575mm**SKE015-4**
Nurse StoolChrome Plated Structure with PU Leather
High Density Sponge Seat
Size: 440×400×450-590mm**SKE013-11**
Lab ChairChrome Plated Structure with Polyurethane Foam Seat
Size: 330×220×450-520mm**SKE013-10**
Lab ChairChrome Plated Structure with Polyurethane Foam Seat
Size: 430×435×575mm**SKE017-4**
Stainless Steel Surgical Stool201 Stainless Steel Structure with 430 Stainless Steel Seat
Size: 4255×555×700mm**SKB1A01**Foldable Stretcher
SKA-1C

Material: Aluminium Alloy Pipe & Canvas Fabric



Technical Parameters

Full Size	2100×500×100mm
Folded Size	580×220×210mm
Net Weight	6kg
Load Bearing	159kg



4-fold

SKB1A02Foldable Stretcher
SKA-1B

Material: Aluminium Alloy Pipe & Canvas Fabric



Technical Parameters

Full Size	2200×530×140mm
Folded Size	1100×180×130mm
Net Weight	6kg
Load Bearing	159kg



2-fold

SKB1A04Foldable Stretcher
SKA-1D

Material: Aluminium Alloy & Canvas Fabric



Technical Parameters

Full Size	1840×510×210mm
Folded Size	920×510×100mm
Net Weight	7.5kg
Load Bearing	159kg




IB 4-fold

IA 2-fold

SKB2B01Scoop Stretcher
SKX-B1

Material: Aluminium alloy




Technical Parameters

Full Size	2200×420×70mm
Folded Size	1660×420×70mm
Net Weight	8.5kg
Load Bearing	159kg



SPINE BOARD



SKB1C10

Electric Stair Stretcher
SKW-B2



Technical Parameters

Full Size	1230×500×1400mm
Folded Size	1020×550×210mm

SKB1C01-1

Foldable Transfer Stair Stretcher
SKW-B2



Technical Parameters

Full Size	850×520×1370mm
Folded Size	1100×540×180mm

SKB1C02-1

Manual Stair Stretcher for Downstairs
SKW-B2



Technical Parameters

Full Size	1250×550×1630mm
Folded Size	1250×550×330mm

Ambulance Stretcher Trolley

Material: Aluminium Alloy / Foam Mattress / PVC Cover

SKB039(A)

SKC-2B1



Technical Parameters

Full Size	1900×550×400mm	Net Weight	25.5kg
Back-rest	70°	Load Bearing	159kg

SKB039(E)

SKC-3F5



Technical Parameters

Full Size	1900×540×800mm	Net Weight	45kg
Lowest Size	1900×540×200mm	Load Bearing	159kg
Back-rest	85°		

SKB039(C)

SKC-3B4



Technical Parameters

Full Size	1900×550×920mm	Net Weight	36.8kg
Lowest Size	1900×550×260mm	Load Bearing	159kg
Back-rest	70°		

SKB039(D)

SKC-3B5



Technical Parameters

Full Size	1950×550×870mm	Net Weight	40kg
Lowest Size	1950×550×230mm	Load Bearing	159kg
Back-rest	70°		

SKB037(A)

Stretcher Trolley

SKM-A

Size: 1960×600×810mm

**SKB037(B)**

Patient Trolley

SKM-C

Size: 1960×600×770mm



70°

Back rest

**SKB037(C)**

Patient Trolley

SKM-G

Size: 1970×760×620-900mm



75°

Back rest




35°

Knee rest



280mm

Height Adjustment
620-900mm**SKB038-2**

Patient Trolley

SKM-E

Size: 1900×700×620-900mm



65°

Back rest



280mm

Height Adjustment
620-900mm**SKB038-4**

Patient Trolley

SKM-G

Size: 1900×700×620-900mm



65°

Back rest



280mm

Height Adjustment
620-900mm**SKB041-3S**

Hydraulic Patient Trolley

SKM-G

Size: 2150×790×635-895mm



60°

Back rest



260mm

Height Adjustment
635-895mm**SKB041-1**

Patient Transportation Trolley


SKM-F

70°±10°

Back rest



300mm

Height Adjustment
540-840mm

Technical Parameters

Length	1900mm
Width	730mm
Min. Height	550mm
Max. Height	850mm

Technical Configuration

✓ Crank System	1 set	✓ Directional Wheel	1 pc	✓ Central Locking Break Pedal	2 pcs
✓ Back Rest Gas Spring	1 pc	✓ IV Pole Previsions	2 pcs	✓ Transfer Mattress	1 pc
✓ ABS Side Rail	1 set	✓ IV Pole	1 set	✓ Directional Pedal	2 pcs
✓ 6" Castors with Central Locking	4 pcs	✓ Oxygen Cylinder Holder	1 pc		

SKB041-2

Patient Transportation Trolley


SKM-F

75°±10°

Back rest



300mm

Height Adjustment
630-930mm

15°

Reverse Trendelenburg



15°

Trendelenburg



Technical Parameters

Length	1930mm
Width	640mm
Min. Height	630mm
Max. Height	930mm

Technical Configuration

✓ Hydraulic Pump	2 pcs	✓ Directional wheel	1 pc	✓ Directional Pedal	2 pcs
✓ Back Rest Gas Spring	1 pc	✓ IV Pole Previsions	2 pcs	✓ Transfer Mattress	1 pc
✓ ABS Side Rail	1 set	✓ IV Pole	1 set		
✓ 8" Castors, Central Locking	4 pcs	✓ Central Locking Break Pedal	2 sets		

SKB041-3
Patient Transportation Trolley
SKM-G



Technical Parameters

Length	2100mm
Width	850mm
Min. Height	630mm
Max. Height	930mm

Technical Parameters

✓ USA Hydraulic Pump	2pcs	✓ Directional Wheel	1pc	✓ Ø74mm Wall Bumper Wheel	4pcs
✓ Crank System	1set	✓ Folding Side Rail	1set	✓ Central Locking Break Pedal	2sets
✓ Back Rest Gas Spring	2pcs	✓ IV Pole	1set	✓ Directional Wheel Release Pedal	2pcs
✓ 8" Castors, Central Locking	4pcs	✓ IV Pole Previsions	4pcs	✓ X-ray Translucent Bed Platform	1set
✓ Emergency Push Handle	1set	✓ Mattress	1pc		

SKB041-6
Patient Transportation Trolley
SKM-F



Technical Parameters

Length	1930mm
Width	660mm
Min. Height	510mm
Max. Height	850mm

Technical Parameters

✓ Crank System	1set	✓ Directional Wheel	1pc	✓ Oxygen Cylinder Holder	1pc
✓ Back Rest Air Spring	1pc	✓ Directional Pedal	2pcs	✓ Central Locking Break Pedal	4pcs
✓ Plastic Side Rail	1set	✓ IV Pole Previsions	2pcs	✓ Transfer Mattress	1pc
✓ 6" Castors, Central Locking	4pcs	✓ IV Pole	1set	✓ Safe Lock for Side Rail	2pcs

SKB041-10
Patient Transportation Trolley
SKM-F



Technical Parameters

Length	1880mm
Width	620mm
Min. Height	560mm
Max. Height	890mm

Technical Parameters

✓ Crank System	1set	✓ Directional Wheel	1pc	✓ Central Locking Break Pedal	4pcs
✓ Back Rest Air Spring	1pc	✓ Directional Pedal	2pcs	✓ Transfer Mattress	1pc
✓ Plastic Side Rail	1set	✓ IV Pole Previsions	2pcs	✓ Safe Lock for Side Rail	2pcs
✓ 6" Castors, Central Locking	4pcs	✓ IV Pole	1set	✓ Emergency Push Handle	1set
✓ USA Hydraulic Pump	2pcs	✓ Back Rest Gas Spring	1pc	记录台	1pc
✓ Wall Bumper Wheel	4pcs	✓ Oxygen Cylinder Holder	1pc		

SKB041
Operation Connecting Trolley
SKK-A



Technical Parameters

Length	3650mm
Width	640mm
Min. Height	640mm
Max. Height	970mm

Technical Configuration

✓ Crank System	2pcs	✓ 6" Castors with Central Locking	8pcs	✓ Directional Wheel Release Pedal	4pcs
✓ Back Rest Air Spring	1pc	✓ IV Pole Previsions	2pcs		
✓ ABS Side Rail	1set	✓ IV Pole	1set		
✓ Directional Wheel	2pcs	✓ Central Locking Break Pedal	4pcs		

SKR-AB00

New All-in-One Computer Cart
Size: 490×780×980-1030mm

**SKR-A01**

Simple All-in-one Computer Cart
Size: 650×490×1260mm

**SKR-B00**

Simple All-in-one Computer Cart
Size: 533×455×750-1150mm

**SKR-IB00**

Simple Laptop Cart
Size: 533×455×750-1150mm

**SKR-J10**

Patient Monitor Cart
Size: 505×630×950mm

**SKR-J21**

Endoscopic System Equipment Cart
Size: 665×630×1650mm

**SKR-R06**

Monitor Cart
Size: 510×1000mm

**SKR-R05**

Tablet Cart
Size: 510×830mm

**SKR-R03**

Laptop Cart
Size: 530×520×850mm

**SKR-R04**

Ward Rounds Cart
Size: 530×520×730mm

**SKR-R02**

Mobile Multi-function Information Search Cart
Size: 530×520×1200mm

**SKR-R01**

Information Nursing Cart
Size: 520×620×1315mm



Registered Model: SKM-A

SKR054-ET

Emergency Trolley
Size: 750×475×920mm

**Technical Configuration**

✓ Defibrillator shelf	1pc	✓ IV pole	1pc	✓ Central Lock	1set
✓ Sliding Side Shelf	1pc	✓ CPR Board	1pc	✓ File Box	1pc
✓ Sharp Container	1pc	✓ Oxygen Bottle Holder	1pc		
✓ Dust Baskets	2pc	✓ Power Outlet	1pc		

Registered Model: SKM-B

SKR051-ET

Emergency Trolley
Size: 750×475×930mm
850×520×960mm

**SKR054-AT**

Anesthesia Trolley
Size: 750×475×930mm
850×520×960mm

**SKR054-MT**

Medicine Trolley
Size: 750×480×920mm
850×520×950mm

**SKR-CT750-1**

Clinical Trolley
Size: 750×475×920mm

**SKR-ET625**

Emergency Trolley
Size: 625×475×930mm

**SKR-AT625-1**

Anesthesia Trolley
Size: 625×475×930mm

**SKR-MT625**

Medicine Trolley
Size: 625×475×930mm

**SKR-CT625-1**

Clinical Trolley
Size: 625×475×920mm

**SKR004**

ABS Trolley
Size: 730×480×900mm

**SKR007**

ABS Trolley
Size: 730×480×900mm

**SKR021-8**

Case History Trolley
Size: 750×410×950mm

**SKR022**

Case History Trolley
Size: 450×410×950mm

**SKR001-2**

ABS Trolley
Size: 539×366×845mm
640×456×845mm
760×520×845mm

**SKR002-3**

ABS Trolley
Size: 539×366×855mm

**SK-MT650**

Medicine Trolley
Size: 650×492×1040mm

**SK-ET750**

Emergency Trolley
Size: 750×530×1050mm

**SK-ET75077A**

Emergency Trolley
Size: 750×475×930mm

**SK-AT75077A2**

Anesthesia Trolley
Size: 750×475×930mm

**SK-CT75077C2**

Treatment Trolley
Size: 750×475×920mm

**SKR-ET850**

Emergency Trolley
Size: 850×520×990mm

**SKR-CT600**

Clinical Trolley
Size: 775×635×910mm

**SKR-MT606**

Clinical Trolley
Size: 775×630×1080mm

**SKR-CT232**

Clinical Trolley
Size: 780×475×1000mm

**SKR-CT601**

Clinical Trolley
Size: 775×635×910mm

**SKR-CT670**

Clinical Trolley
Size: 670×450×825-1000mm

**SKR-CT670-1**

Clinical Trolley
Size: 670×450×920-950mm

**SKR-AT670**

Anesthesia Trolley
Size: 670×450×880-960mm

**SKR-ET670**

Clinical Trolley
Size: 670×450×880-960mm

**SKH040**

Laundry Trolley
Size: 800×600×840mm

**SKH040-1G**

Laundry Trolley
Size: 800×600×845mm

**SKH027**

Laundry Trolley
Size: 900×450×1050mm

**SKR-WC561**

Waste Collecting Trolley
Size: 800×590×890mm



SKH006C

Instrument Trolley
[Disassembly structure]
Size: 670×450×900mm

**SKH006**

Instrument Trolley
Size: 670×450×900mm

**SKH001**

Treatment Trolley
Size: 670×450×900mm

**SKH002**

Dressing Trolley
Size: 670×450×900mm

**SKH006-1**

Treatment Trolley
Size: 670×450×900mm

**SKH004**


Instrument Trolley
Size: 670×450×900mm

**SKH006-5**

Treatment Trolley
Size: 460×460×860mm

**SKH005**

Fan-shaped Instrument Trolley
Size L: 1550×460×900mm

**SKH029**

Dressing Trolley
Size: 670×450×900mm

**SKH018**

Anesthesia Trolley
Size: 600×400×850mm

**SKH032**

Infusion Treatment Trolley
Size: 670×450×1700mm

**SKH012**

Insulated Food Cart
Heat store type
Size: 1030×550×980mm

**SKH021**

Emergency Trolley
Size: 700×400×850mm

**SKH019**

Emergency Trolley
Size: 700×400×850mm

**SKH028**

Medicine Trolley
Size: 670×450×1000mm

**SKH007-1**

Stainless Steel Trolley
Size L: 1000×500×900mm
H: 880×430×860mm
S: 760×360×820mm

**SKH013**

Stainless Steel Trolley
for Medical Record Holder (40 shelves)
Size: 630×430×900mm

**SKH014**

Stainless Steel Trolley
for Medical Record Holder (20 shelves)
Size: 340×400×900mm

**SKH070**

Work Table
Size: 1200×600×800mm

**SKH077**

Work Table
Size: 1200×800/1000×1600mm

**SKH078**

Shelf With Four Layers
Size: 1200×600×1800mm

**SKH105**

Wet Linen Trolley
Size: 815×610×750mm

**SKH048**


Oxygen Bottle Trolley
Size: 1435×360×1040mm

**SKH088**


Swab Stand
Size: 900×500×1500mm

**SKH036**


Stainless Steel Washing Sink
Size: 1200(W)×600(D)×1100(H)mm



Size: 600(W)×600(D)
×1100(H)mm



Size: 1800(W)×600(D)
×1100(H)mm



Size: 2400(W)×600(D)
×1100(H)mm

**SKH041 (5)**

Wash Basin
Size: Φ520×810-1300mm



SKH041 (6)

SKH038-2

Mayo Trolley
Size: 500×350×870-1300mm



SKN001

Stainless Steel Tank
 $\Phi 170\text{mm}$ $\Phi 200\text{mm}$ $\Phi 230\text{mm}$
 $\Phi 260\text{mm}$ $\Phi 300\text{mm}$

**SKN002**

Stainless Steel Cream Gauze Can
 $\Phi 80 \times 80\text{mm}$ $\Phi 90 \times 90\text{mm}$
 $\Phi 100 \times 100\text{mm}$ $\Phi 120 \times 120\text{mm}$
 $\Phi 141 \times 125\text{mm}$ $\Phi 160 \times 154\text{mm}$

**SKN011**

Stainless Steel Sterilization Tray
 Without holes: With holes:
 $150 \times 100 \times 45\text{mm}$ $150 \times 100 \times 45\text{mm}$
 $200 \times 110 \times 45\text{mm}$ $200 \times 110 \times 45\text{mm}$
 $250 \times 160 \times 45\text{mm}$ $250 \times 160 \times 45\text{mm}$
 $300 \times 200 \times 50\text{mm}$ $300 \times 200 \times 50\text{mm}$

**SKN012**

Stainless Steel Tray
 $300 \times 200 \times 48\text{mm}$ $350 \times 250 \times 48\text{mm}$
 $400 \times 300 \times 48\text{mm}$ $450 \times 300 \times 48\text{mm}$
 $350 \times 500 \times 48\text{mm}$

**SKN021**

Stainless Steel Needle Sterilizer
 Size: 77holes 105holes

**SKN022**

Stainless Steel Medicine-change Bowl
 $\Phi 180\text{mm}$ $\Phi 200\text{mm}$ $\Phi 220\text{mm}$
 $\Phi 240\text{mm}$ $\Phi 260\text{mm}$ $\Phi 280\text{mm}$

**SKN031**

Stainless Steel Kidney Dish
 Shallow: Deep:
 $L204 \times 130\text{mm}$ $L243 \times 125\text{mm}$
 $M190 \times 116\text{mm}$ $M203 \times 102\text{mm}$
 $S178 \times 104\text{mm}$ $S165 \times 80\text{mm}$

**SKN032**

Stainless Steel Sterilized Bottle/
 Medicine Cup
 Sterilized bottles:
 $\Phi 50 \times 103\text{mm}$ $\Phi 50 \times 93\text{mm}$ $\Phi 50 \times 83\text{mm}$ $\Phi 50 \times 53\text{mm}$
 Medicine cups: $\Phi 44 \times 40\text{mm}$ $\Phi 50 \times 50\text{mm}$

**SKN041**

Stainless Steel Soak Forceps Bucket/
 Measuring Cup/Enema Bucket
 Measuring cup: Enema bucket:
 1L $\Phi 100 \times 150\text{mm}$ 1L $\Phi 100 \times 150\text{mm}$
 2L $\Phi 130 \times 190\text{mm}$ 2L $\Phi 130 \times 190\text{mm}$

**SKN042**

Stainless Steel Forceps Bucket
 Inverted cone: Straight:
 $\Phi 80 \times 198\text{mm}$ $\Phi 80 \times 200\text{mm}$
 $\Phi 64 \times 158\text{mm}$ $\Phi 64 \times 160\text{mm}$
 $\Phi 45 \times 113\text{mm}$ $\Phi 45 \times 113\text{mm}$

**SKH090S**

Cabinet
 Size: $900 \times 400 \times 750\text{mm}$

**SKH090M**

Cabinet
 Size: $900 \times 400 \times 1050\text{mm}$

SKH090L

Cabinet
 Size: $900 \times 400 \times 1800\text{mm}$

**SKH092**

Cabinet
 Size: $900 \times 400 \times 1800\text{mm}$

**SKH093**

Cabinet
 Size: $900 \times 400 \times 1800\text{mm}$

**SKH098-2**

Locker
 Size: $900 \times 500 \times 1800\text{mm}$

**SKH098-3**

Locker
 Size: $900 \times 500 \times 1800\text{mm}$

**SKH098-4**

Locker
 Size: $900 \times 500 \times 1800\text{mm}$

**SKH098-6**


Locker
 Size: $900 \times 500 \times 1800\text{mm}$

**SKH098-9**

Locker
 Size: $900 \times 500 \times 1800\text{mm}$

**SKH059**

Single-sided
 Medicine Shelf
 Size: $800 \times 460 \times 2000\text{mm}$

**SKH060**

Rotary
 Medicine Shelf
 Size: $920 \times 1780\text{mm}$

**SKH061**

Adjustable Component
 Medicine Shelf
 Size: $1025 \times 500 \times 2000\text{mm}$

**SKH062**

Two-sided Adjustable
 Medicine Tray
 Size: $800 \times 860 \times 2000\text{mm}$



SKH059-1
 Size: $800 \times 660 \times 2000\text{mm}$

NOTE

All the images, specifications and illustrations belong to SAIKANG's proprietary products. Any form of imitation or copy is a violation of the rights of our company. The actual color of products maybe have color difference as shown in catalogue. Products appearance maybe have difference if when technique upgraded.
Copyright © 2022 SAIKANG All Rights Reserved.

全国服务热线: 400 828 7009
www.jssaikang.com

赛康 (美国)
Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, Ca 91748
Tel: +1 949648-2776

赛康 (东非)
Saikang Medical Solutions (East Africa) Ltd.
Vinej Complex Unit 3 Mombasa Road,
PO Box 14262-00800, Nairobi
TEL: +254 733-704 690

赛康 (总部)
江苏赛康医疗设备股份有限公司
江苏省张家港市乐红路35号

赛康 (成都)
成都赛康医疗设备有限公司
成都市武侯区西部智谷B区5栋1单元5F

赛康 (英国)
Uk Healthward International Limited
Owner(s) Address DEPT 111,
196 High Road Wood Green,
London, N22 8HF,
United Kingdom
TEL: +44 2032399738

赛康 (中东)
Saikang Medical Supplies (Mid East) LLC
XQ22-S03 England Cluster,
International City, Dubai UAE
TEL: +20 128 8962497



微信公众号



赛康官网

本目录中的内容 (图像、技术说明、插图) 是赛康医疗的专有产品宣传资料, 任何形式的仿制和任何媒介的复制, 都是对该公司权利的侵害。由于印刷误差, 图片与实物产品的颜色存在差异, 以实物为准。产品会随设计更新、改良、技术参数以及外观的变化而变化。敬请谅解。图样中所采用的人物照片为其肖像权购买, 如侵犯肖像权, 请与设计方联系, 并提供真实有效的所有权证明, 以便维护双方权益。

Copyright © 2017/2022 SAIKANG All Rights Reserved

SAIKANG
医院家具

www.jssaikang.com

C21/22-1

SAIKANG[®] MEDICAL FURNITURE

医院家具

江苏赛康医疗设备股份有限公司
www.jssaikang.com

股票代码: 870098

SAIKANG[®] Medical



SAIKANG

与时间赛跑·助生命健康

RACE AGAINST TIME, HEALTHIER YOUR LIFE

www.jsaikang.com



RACE AGAINST TIME HEALTHIER YOUR LIFE

与时间赛跑·助生命健康

www.jsaikang.com

CONTENTS 目录

Saikang Medical	门诊部	住院部	行政中心	医技部	Saikang Service
01 实力赛康	07 导医台	35 护士站	63 密集柜	73 检验科	77 服务网络
03 科技赛康	09 收费挂号	39 治疗室	65 院长办公区	75 超声科	
	11 自动化药房	41 处置室	66 主任办公区		
	13 西药区	42 换药室	67 医生办公区		
	15 中药区	43 配液室	68 会议室		
	17 诊室	44 污洗间	71 多功能厅		
	27 诊区配套	45 病房			
	31 输液室	57 医护值班/休息室			
	33 雾化区	59 医护更衣室			
	34 洗眼中心				

STRENGTH OF SAIKANG

实力赛康

20

20年专注医疗护理设备设计、生产、服务及技术咨询于一体的科技型企业。

150

地处以上海为经济中心的150公里经济辐射圈之苏州-张家港。

10000000

赛康累计向全球输送100多万套产品，产品遍布亚洲、中东、北美和非洲等地区。

400

全球雇员400余人，高中级职称50余人，技师100余人。

189

拥有完善的质量体系和专业的检测团队，以及各类自动化生产设备189台。每个生产环节和生产过程都有严格的标准控制，确保每个零件和产品的质量，以及可追溯的标准管理。

50000

一期占地面积30,000多平方米，二期20,000多平方米现代化生产基地3年内投产。

18

赛康提供完善的医疗解决方案，拥有18大品类。

562

安全、高效、易用的整体解决方案，包括急救、护理和ICU病房等500多项系列产品。

29

先后通过德国TUV的欧盟CE安全认证、ISO9001/13485质量体系认证、美国FDA市场许可等29项认证。

131

累计申请各项专利131项，不断领先定义产品质量和服务效率的新标准。

149

公司成立以来，产品与解决方案已应用于全球超过149个国家及地区。

1000

赛康产品致力于打造现代化的护理环境，累计服务超1000家医疗与护理机构。



LEADING BRAND OF MEDICAL FURNITURE

整体医护设备一站式解决方案

医疗康复设施中的病人常常不得不面对巨大的挑战。许多人遭受了疾病或事故，面临着重新训练自己的身心，以及必须适应新环境。在这种情况下，积极的气氛可能有着决定性的帮助。床和家具应该为病人和护理人员提供最好的支持。赛康重视医疗护理解决方案，为医疗机构及康复中心帮助病人实现快速和长期恢复的完美选择。

赛康创始于2002年，是全球领先的医疗设备与解决方案供应商。赛康总部设在中国苏州，在中国17个省市自治区均设有办事处，全球雇员超过400名，形成了庞大的全球研发、营销和服务网络。赛康融合创新，紧贴临床需求，帮助世界各地人们改善医疗条件，降低医疗成本。目前，赛康的产品与解决方案已应用于全球超过149个国家及地区，中国超千家医疗机构和37%以上的三甲医院。

TECHNOLOGY OF SAIKANG

科技赛康

赛康坚持提供高品质的产品。在欧美各大洲市场新推出产品，对“中国制造”产品的信心。赛康先后通过美国FDA、CE、ISO 9001/13485质量管理体系认证、欧盟CE安全认证、美国FDA（美国食品药品监督管理局）市场准入许可。

Alibaba.com Assessed Supplier

Jiangsu Saikang Medical Equipment Co., Ltd.
江苏赛康医疗设备有限公司

中国制造网 认证供应商
AUDITED SUPPLIER
Jiangsu Saikang Medical Equipment Co., Ltd.
江苏赛康医疗设备有限公司

Adams China.com
中国医药网

苏州名牌产品
江苏省名牌产品



Outpatient

Inpatient

Administration Center

Medical Technology

SAIKANG[®]

MEDICAL FURNITURE

医院家具一站式智造商

满足医院的使用需求、安全耐用、操作简单、易于清洁、抗感染，能创造良好的使用环境，较高的使用寿命。但作为医用家具，能否满足人文需求尤显重要。人文需求满足了能减轻病人的心理压力，起到安抚病人情绪和舒缓心情的作用。

安全性

符合院感管理，洁净防污，材质环保，施工过程环保，结构稳固安全

耐久性

材料品质保障，结构坚固耐用

人性化

功能人性，提升舒适性，提升医院服务和运作效率，合理规划设计与智能化应用

经济化

材料、设计、生产、施工，适应医院成本控制，践行节俭社会价值观

智能化

更加便捷、舒适，未来医疗家具的发展趋势

GUIDANCE TABLE

导医台



台面

优质复合亚克力医用人造石

框架

20mm厚一级电解板框架, 1.0mm厚一级电解板饰面, 进口抑菌粉末静电喷涂

脚线

304不锈钢

五金

阻尼导轨, 不锈钢缓冲门铰, 三折一体式拉手

材质: 面板: 亚克力板 柜体: 金属板



TOLL REGISTRATION

收费挂号

模块化屏风隔断

主体采用国际1.0mm医用电解板，经数控剪板、冲压、折弯、激光切割、焊接等一系列程序一体成型，外型设计前卫美观，流线层次分明。

抑菌粉末喷涂，抗大肠杆菌、金黄色葡萄菌等多种菌种。

台面采用12.7mm人造石。



饰面材质

全钢制饰面材料经静电喷涂处理

内部结构

全钢制框架，玻璃四边收口，防震处理

工艺处理

工业化4.0设备加工，圆角收边收口工艺

施工周期

现场安装2个工作日内完成+运输周期

统一性

工厂机械化、标准化生产，确保医院形象统一

快捷性

较装修施工周期缩短一半以上，提高整体施工效率

安全性

全钢制材料结构，确保安全性能达到要求

灵活性

独特组装方法，现场预留尺寸，不影响其他施工环节

配套性

侧板、台架、灯箱等全部自行生产，配套性能完整

环保性

大部分为金属可回收材料，现场施工无污染

重组性

全部为框架结构，现场组装，重复使用率100%



AUTOMATED PHARMACY

自动化药房

全流程体系的自动化解决方案，从药品的一级缓存库到药品的调剂、传送、前台配发均有相关自动化设备处理。



自动化二级缓存区

采用机械手存取系统，将药品自动存储在自动化二级库中，根据药品消耗情况，自动生成请领计划并自动发送。



高速发药机

用于高频次大批量药品存储与发放的专用发药系统。针对此类药品，采用高效的存储方式，单槽位每秒3盒的发药速度，解决以往设备屏蔽大处方的瓶颈，有效提高了药房的自动化率。



自动发药机

纯盒装处方药品直接通过自动传输系统发送至前台，配合前台发药支持系统，发药人员无需走动。



混合处方区域

用于盒装药品和异型包装药品混合处方的调剂，调剂人员无需大范围走动，将药品集合核对后，通过整处方传递系统传送至发药前台。



WESTERN MEDICINE

西药区

柜体

电解钢板，钢板厚度 $\geq 1.0\text{mm}$ ；主体钢结构，采用 2.0mm 电解钢板和 2.0mm 钢管结合构造，柜门与抽屉采用隐藏式拉手，折弯一体成型

台面

复合亚克力树脂 $\geq 28\%$ +铝粉，易清洁、耐高温、耐冲击、抗污垢

导轨

采用专业级医用导轨，具有顺畅、静音、防倾覆功能

表面工艺

采用静电环氧树脂粉末喷涂


踢脚线

采用304不锈钢材质



SKH059

单面药架
800×460×2000mm



SKH059-1

双面药架
800×860×2000mm



SKH062

双面可调式药架
800×860×2000mm



SKH060

旋转式药架
4820×1790mm



SKH061

组合式药架
1025×500×2000mm



CHINESE MEDICINE

中药区



SK-YG200

实木中药柜
4025x640x2620mm



SK-YG201

实木中药柜
2307x740x2300mm




SK-YG202

实木中药柜
2278x725x2115mm



SK-YG001

单面中药柜



SK-YG002

双面中药柜



SK-YG003

单面中药调剂柜



SK-YG004

双面中药调剂柜



ORDINARY CLINIC

普通诊室

运用整体医院家具设计风格，以简洁耐用的家具材料为主，局部可采用高档家具作为点缀色彩，色调柔和以平和病人的情绪，营造亲切宜人并相对私密的诊疗环境。



PEDIATRIC CLINIC

儿科诊室

运用整体医院家具设计风格，发挥色彩的辅助治疗作用。主要针对病婴幼儿，不同年龄的儿童特征，营造活泼可爱、舒适温馨、安全现代的诊疗场景氛围。



GYNECOLOGY CLINIC

妇科诊室

运用整体医院家具设计风格，色彩方面以暖色为主，局部采用粉色调进行点缀，营造和谐平静的温馨氛围。设备上考虑手术诊断检查等多种医疗功能使用。



DENTAL CLINIC

口腔科门诊

口腔诊室是口腔专科工作场所，医生一般需要借助专业的牙椅对患者进行检查和治疗，口腔诊室可以是单独房间，也可以是开放空间分出隔断诊室。



CHINESE MEDICINE CLINIC

中医门诊

设计采用新中式风格，设计配套诊室检查、隔断、问诊桌椅等，同时兼容医疗诊断设备。



SUPPORTING CLINIC

诊区配套




SKE008-1



不锈钢款



SKE008

SKE006-1
1655×675×870mm(L×W×H)SKE006-2
1765×675×905mm(L×W×H)

SUPPORTING CLINIC

诊区配套


根据客户定制



INFUSION ROOM

输液室

根据患者病情，设输液椅、输液床供患者使用。选材上采用环保、坚固耐用、耐褪色、抗酸碱腐蚀性的材料，能使用消毒液清洗。金属表面光滑，避免划伤患者。



NEBULIZER THERAPY AREA

雾化区

主要用于急诊、儿童患者雾化的治疗。采用整体雾化室家具解决方案，平面布局严格按照国家雾化区感控标准。家具设计及色彩的运用充分考虑恐惧、烦躁的心理特征以及病人隐私心理需求。



BABY WASH CENTER

洗婴中心

新生儿科对医疗环境安全、人为关怀等方面要求特殊，注重人性化设计。所配套的家具达到无毒害舒适安全，色彩欢快具有趣味性的功能。



NURSE STATION

护士站

医院服务台是体现医院或者医疗机构形象和档次的门户所在，设计要充分考虑达到亲切优雅、拉近距离、融洽和谐的概念。



NURSE STATION

护士站



TREATMENT ROOM

治疗室

治疗室一般靠近护士站布局，作为准备输液、配药备车、准备辅助器械等的科室。上柜及中柜空间延伸，让存储空间最大化利用，整体的动线流程布局，保证了医护工作的快、准、精。



DISPOSAL ROOM

处置室

处置室位于治疗室的近旁，作为医疗废物的收纳房间，对医疗垃圾、生活垃圾、可回收垃圾做科学的分类处理。设有垃圾柜，治疗柜，以及放置各种医疗器具的操作台。



DRESSING ROOM

换药室

换药室的空间设计应布局合理分区明确无菌物品和污染物品分开放置。



DOSING ROOM

配液室

配液室设计在人员流动较少的安静区域，便于医护人员沟通和成品的运送，设置地点远离各种污染源。配液室的设计，充分考虑到洁净区和辅助工作区。



SLUICE ROOM

污洗间

污洗室设计靠近污物出口处，可设计衣物清洗设施和投放开口，还可配置清洁车、置物架、垃圾桶等。



WARD

普通病房

普通病房的设计运用整体医院家具设计风格。考虑提供病人休息康复治疗的设备。为了便于清洁及护理行为的需要，病床三面临空。每一个床单元应配备床头柜、陪护椅、移动餐桌、隔帘、设备带等基本设施。



VIP ROOM

VIP病房

VIP房间的设计，充分考虑到使病房的配置更加人性化，适当加入沙发、茶几、家用电器等家用设施，能够提供温馨的康复治疗环境。



ICU

ICU病房

ICU通常采用开放式多床布置方式，但从感染控制、防止交叉感染、保护患者与医护人员安全的角度出发，应鼓励在人力资源充足的条件下，尽量多设计单间或分隔式病房，尤其是以内科患者为主的ICU。



ORTHOPEDICS

骨科病房

运用整体医院家具设计风格，病床的设计上充分考虑到了骨科患者的康复，需要配有骨科支架、牵引架等医疗设施。



GYNECOLOGICAL

妇科病房

运用整体医院家具设计风格，色彩上可以考虑淡紫色和淡粉色为主色调，舒缓紧张的情绪，让人感觉优雅而宁静。



PEDIATRIC

儿科病房

儿科对于医疗环境安全、人为关怀等方面要求环境特殊，注重人性化医疗环境配备，必须考虑到卡通化的优化环境，安全防护方面要有防坠床设施以及防撞措施，配套的家具需达到无毒害、舒适、安全、色彩欢快等特点。



MEDICAL DUTY/LOUNGE

医护值班/休息室

值班室也称休息室是医护人员公休期间临时休息的场所。设置在病区相对安静的地方靠近更衣室。既能保障医护人员的休息，又能在其他工作人员呼叫的同时，比较短的路径，最快的到达病房。



CHANGING ROOM

医护更衣室



更衣区临床医护人员每天要接触大量病患。长期工作在细菌病毒密集的环境中，医护的卫生行为是医院感染控制中非常重要的一个环节。更衣卫浴靠近值班室的合理设置，使得功能集中，方便医护人员使用节省空间，互不干扰。





SKH057A
18门彩喷鞋柜
900×400×1800mm



SKH058-5A
彩喷更衣柜
900×400×1800mm



SKH051A
彩喷器械柜
900×400×1800mm



SKH052A
彩喷针剂柜
900×400×1800mm



SKH055A
上柜操作台
1600×500/300×1750mm



SKH040-12
三联脚踏污物柜



SKH040-11
存车柜



SKH036-15
顶柜/水池柜

SKH090S
文件柜
900×450×750mm



SKH090M
文件柜
900×450×1050mm



SKH090L
文件柜
900×450×1800mm




SKH098-2
更衣柜
900×500×1800mm



SKH098-3
更衣柜
900×500×1800mm



SKH098-4
更衣柜
900×500×1800mm



SKH098-6
更衣柜
900×500×1800mm



SKH098-9
更衣柜
900×500×1800mm



SKH091
文件柜
900×450×1800mm



SKH092
文件柜
900×450×1800mm




SKH093
文件柜
900×450×1800mm



MOBILE FILING SHELF

密集柜





DIRECTOR OFFICE

主任办公区

SKZ311

办公桌

桌子: 1510x810x750mm
 横柜: 1800x550x470mm



SKZ313

办公桌

桌子: 2200x900x750mm
 横柜: 1300x650x400mm



SKZ312

办公桌

1200x600x750mm



SKZ314

办公桌

1800x1200x1040mm



DEAN OFFICE

院长办公区

SKZ432

办公桌

桌子: 2400x950x750mm
 横柜: 2000x800x445mm



SKZ433

办公桌

1800x800x750mm



SKZ434

办公桌

桌子: 1800x800x750mm
 横柜: 1800x600x445mm





DOCTOR OFFICE

医生办公区

医生办公室是医生办公学习交流研讨的场所，设计考虑在大空间内用矮隔断划分空间，一般设置小型会议桌，方便满足医生研讨之用，同时也可供实习医生办公学习。功能方面可充分考虑办公区、交流区、储物区等医院家具。

SKZ107

办公桌

1400×600×750mm



SKZ110

办公桌

1200×700×750mm



SKZ111

办公桌

1400×700×750mm



SKZ191

办公桌

1200×600×740mm



MEETING ROOM

会议室

设计充分考虑满足会议示教、远程会诊等要求，可提供会议区、演讲区、辅助区等家具模块。

SKZ701

会议桌



SKZ702

会议桌



SKZ801

会议椅





SKE702
医师椅

480×430×460-560mm



SKE568
医师椅

630×550×460-560mm



SKE714
医师椅

690×620×470-550mm



SKE705
医师椅

490×430×440-540mm



SKE711
医师椅

500×490×500-600mm



SKE703
医师椅

630×530×440-540mm



SKE054-2
医师椅

570×520×410-620mm



SKE055-8
医师椅

480×470×900mm



SKE806
折叠椅

折叠尺寸: 190×910mm
展开尺寸: 525×600×795mm



SKE713
医师椅

610×590×950mm



SKE052
会议椅

480×470×780mm



SKE063
医师椅

580×450×980mm



SKE015-2
护士椅

φ340×490-660mm



SKE016
护士椅

φ345×485-800mm



SKE013
护士椅

φ330×445-575mm



SKE013-2
护士椅

430×390×430-570mm



SKE013-11
护士椅

330×230×450-580mm



SKE017-4
护士椅

φ295×500-700mm



SKE013-10
护士椅

φ350×435-575mm



SKE802
牙科椅

φ350×450×450-560mm



SKE020-1
脚踏凳

490×410×170/340mm



SKE020-10
脚踏凳

400×300×90mm



SKE020-4D

405×200×220mm

SKE019
小方凳

380×280×475mm



SKE019-2
小方凳

340×270×440mm



FUNCTION ROOM

多功能厅



SKE048
礼堂椅



SKE164
礼堂椅



SKE161
礼堂椅



SKE162
礼堂椅



SKE165
礼堂椅



SKE166
礼堂椅



SKE170
礼堂椅



SKE171
礼堂椅



LABORATORY

检验科

对于生物、医药实验室、检验科而言，家具的清洁和消毒是首先要考虑的功能，要求其家具具有易清洁、易消毒，并能够较长时间保持洁净的特点。



台面

12.7mm实芯理化板，边缘加厚至25.4mm

柜架

优质C型矩形钢管制作，表面静电粉末喷涂，耐酸碱、耐腐蚀

柜体

采用15mm环保E1级中纤板，2.0mmPVC封边

五金配件

一字暗装拉手、实验室专用铰链、三节承重静音滑轨

PP水盆

三叉水龙头、PVC给排水管、试剂架费用另计



ULTRASOUND

超声科





SERVICE NETWORK

服务网络

稳定的产品质量和诚信可靠的供货信誉，赢得了国内外知名企业的青睐。公司产品遍布欧美、东南亚等地区，面向全球。伴随现代信息化的销售网络，使公司大步迈向专业化、规模化发展的快车道。公司期待与您携手共进，共创美好未来。



赛康
累计向全球输送超
1000000
套产品



7×24
小时服务热线
客户问题及时响应

24hrs: 0086-18021231901
www.saikangmedical.com

SAIKANG (Mid East Branch)
Saikang Medical Supplies (Mid East) LLC
Muscat
PO Box 53200, Abu Dhabi, UAE
Tel: +971 585710652

SAIKANG (Russia Branch)
LLC "SAIKANG"
454004, Russian Federation, Chelyabinsk,
Universitetskaya Naberezhnaya
TEL: +7 (917) 303-00-06

SAIKANG (South America Branch)
Saikang Medical Technology, Inc
W3400 Corrientes, Argentina
TEL: +55 (11) 2758-0438
Rua Universal, Hollywood,
São Bernardo do Campo, Brazil
TEL: +55 (11) 2758-0438

SAIKANG (East Africa Branch)
Saikang Medical Technology Ltd
Mezzanine Floor, Taj Tower,
Upper Hill, Nairobi, Kenya
TEL: +254714084804

SAIKANG (South East Asia Branch)
Saikang Medical Technology Ltd
Sala Thammasop Subdistrict, Thaw
Wattana District
Bangkok, Thailand
TEL: +66-642741046

SAIKANG (China Headquarter)
No.15 Lehong Road, Zhonggang Cite
Jiangsu Province, China
TEL: +86-512-58525850
E-mail: inquiry@saikang.com

SAIKANG (UK Branch)
Uk Healthward International Limited
DEPT 111,
196 High Road Wood Green,
London, N22 8HH,
United Kingdom
TEL: +31 (0) 202111106



Factory Tour

Website



SAIKANG OPERATING ROOM

022/23-1

SAIKANG[®] OPERATING ROOM

JIANGSU SAIKANG MEDICAL EQUIPMENT CO., LTD.
www.saikangmedical.com

Stock Code: 870098





SAIKANG
LEADING BRAND OF CHINA
MEDICAL FURNITURE



LEADING BRAND OF CHINA
MEDICAL FURNITURE

RACE AGAINST TIME
HEALTHIER YOUR LIFE



CONTENTS

01	Strength of Saikang	07	Operating Table	55	Service Network
03	Technology of Saikang	19	Obstetric Bed		
		27	LED Operation Lamp		
		43	Medical Bridge, Ceiling Pendant		
		51	Accessories		

STRENGTH OF SAIKANG

20

20 years focusing on design, production, service and technical consultation of medical care equipment.

150

Located in Zhangjiagang City, 150km nearby Shanghai Economic Center.

1000000

Saikang delivered 1,000,000 sets of products globally. Products all over Asia, Middle East, North America and Africa.

400

There are more than 400 employees worldwide, more than 50 qualified engineers, and more than 100 trained technicians.

189

Equipped with 189 sets of various automatic production facilities. Saikang follow strict standard control system on every production process, enable the traceable of each part.

50000

The finished phase I production area covers more than 30,000 square meters, and phase II will be more than 20,000 square meters of advanced production area.

18

Saikang offer comprehensive medical solutions with 18 product categories.

562

Safe, efficient and easy-to-use medical solution including emergency nursing and ICU wards of more than 500 series of products.

29

With CE certification, ISO9001/13485 quality system certification, U.S. FDA and other 29 certifications.

131

Ownership of 131 patents, continuously leading new standards in defining product quality and service efficiency.

149

Since the establishment of the company, products and solutions have been used in more than 149 countries.

1000

Saikang products are committed to creating a modern nursing environment, accumulatively serving more than 1,000 medical and health care organizations.



ONE-STOP SOLUTION FOR OVERALL MEDICAL EQUIPMENT

Patients in medical rehabilitation often have to face huge challenges. Many people who suffering from diseases or accidents are faced with retraining their minds and bodies and must adapt to the new environment. In this case, the positive atmosphere would offer decisive help. Hospital beds and furniture should provide the best support for patients and nurses. Considering this case, Saikang devoted to provide medical care solutions and perfect choices for medical institutions and rehabilitation centers, to help patients achieve rapid or long-term recovery.

Saikang was founded in 2002, now is a leading global provider of medical furniture and solutions. Head office locates in Suzhou, and branches in US, UK, Dubai, Morocco and Kenya, having 400 employees worldwide. Saikang has formed a large global network of research and development, marketing and service. Currently, Saikang products and medical solutions have been used in more than 108 countries and regions in the world.

TECHNOLOGY OF SAIKANG

SAIKANG insists on providing high quality products. All products have been "Made in China" products in accordance with the ISO9001/13485 quality system certification, and all products conform to EU safety certification and CE mark, and have been approved for market access license.

Alibaba.com Assessed Supplier

Jiangsu Saikang Medical Equipment Co., Ltd.
江苏赛康医疗设备有限公司

中国制造网 认证供应商
AUDITED SUPPLIER
Presented to
Jiangsu Saikang Medical Equipment Co., Ltd.
江苏赛康医疗设备有限公司

Made in China
中国 制造





SAIKANG OPERATING ROOM



A501
Electric-Hydraulic
Operating Table
SKT-B

A501 is an advanced, multi-functional electric hydraulic operating table. Its table top is made of high-strength synthetic board and can transmit X-rays. In addition, the platform can be equipped with a 180° rotation function to realize head and leg interchange, which is applicable for C-arm perspective and placement of some special positions.

15° Reverse Trendelenburg	0° Trendelenburg	15° Lateral Tilt left/right	40°/90° Head Plate up/down	85°/40° Back Plate up/down	15°/90° Leg Plate up/down
80mm Kidney Bridge Up	360mm Horizontal Sliding C-Loop Arm	180° Horizontal Rotation	95° Flex	220° Reflex	



Technical Parameters

External Length	2125mm
External Width	550mm
Table Height Min.	600mm
Table Height Max.	950mm

Technical Configuration

✓ Bed Frame	1pc	✓ Arm Holder with Cushion	2pcs	✓ Kidney Bridge Handle	1pc
✓ Mattress	1set	✓ Leg Support with Cushion	2pcs		
✓ Anesthesia Screen Frame	1pc	✓ Foot Support with Cushion	2pcs		
✓ Shoulder Holder with Cushion	2pcs	✓ Hand Controller	1pc		



A106-2

Electric-Hydraulic Operating Table SKT-B

A106-2 is a highly automated electric-hydraulic operating table, which can meet the needs of various surgical operations. The unique design of a super-low position of 500mm, aimed for the brain surgery and ophthalmic surgery.



Technical Parameters

External Length	2080mm
External Width	550mm
Table Height Min.	500mm
Table Height Max.	800mm

Technical Configuration

✓ Anaesthesia Screen	1pc	✓ Foot Support	1pc	✓ Fixing Clamp	7pcs
✓ Hand Remote	1pc	✓ Mattress	1set	✓ Power Line	1pc
✓ Kidney Bridge Handle	1pc	✓ Shoulder Support	1pair		
✓ Arm Support	1pc	✓ Long Fixing Clamp	1pair		



A302

Electric Operating Table SKL-A

A302 is a highly automated electric Multi-functional operating table. Its elevation, lowering, Trendelenburg, reverse Trendelenburg, lateral tilt, back folding, vertical translation and leg folding are all electrically operated. The operating table is equipped with a digital display hand controller to make operation easy and intuitive. With imported large diameter castors, smooth movement, stable brake system. The standard dual control system guarantees the safety of operation.



Material

The base material for the table is high thickness 304 stainless steel, seamless welding by robot, with modern appearance and more stable.



Mattress

Polyurethane foam mattress, waterproof and seamless, avoid liquid in. The connection between table and mattress is buckle joint, easy to disassemble and install, easy to clean.



Head Plate

Foldable head plate, support for more type of surgical operations.



Leg Plate

The up and down folding actions of the leg plate are all completed by electric. The leg plate is removable, effortless and smooth.



Hand Controller

All functions controlled by both hand controller or control panel on table base column. With position LCD display on hand controller, easy to operate.



Battery Could Support

Battery could support about 80-120 complete operation cycles. Working time could be last for one week according to usage and help to keep the operating table running in the event of power failure.

Technical Parameters

External Length	2150mm
External Width	555mm
Table Height Min.	700mm
Table Height Max.	1000mm

Technical Configuration

✓ Bed Frame	1 pc	✓ Arm Holder	2 pcs	✓ Sliding Block	8 pcs
✓ Mattress	1 set	✓ Leg Support	2 pcs	✓ Spanner	5 pcs
✓ Anesthesia Screen Frame	1 pc	✓ Foot Support	2 pcs		
✓ Shoulder Holder	2 pcs	✓ Hand Controller	1 pc		

A301

Electric Operating Table
SKL-A

A301 is a multi-functional electric operating table. Its mattress is made of polyurethane molded material, can be sterilized and cleaned at high temperature, and there are no gaps so to avoid bacterial growth. The platform is made of high-purity phenolic resin board applicable for C-arm, and can be slid up to 350mm in one direction. The X-ray photography can carry out the whole body radiography examination without moving the patient, no blind spots or shadows.



Technical Parameters

External Length	2150mm
External Width	555mm
Table Height Min.	750mm
Table Height Max.	1000mm

Technical Configuration

✓ Bed Frame	1pc	✓ Arm Holder	2pcs	✓ Sliding Block	8pcs
✓ Mattress	1set	✓ Leg Support	2pcs	✓ Spanner	5pcs
✓ Anesthesia Screen Frame	1pc	✓ Foot Support	2pcs		
✓ Shoulder Holder	2pcs	✓ Hand Controller	1pc		

A100-4

Electric Operation Table
SKL-A

A100-4 is an electric operating table that can meet the needs of various surgeries. It is driven by low-voltage DC motor, with safe and stable performance. Super wide table top, super long sliding distance. It can be applicable for C-arm. It has a high level of automation, low noise, strong reliability, and imported original parts for key parts. It is a perfect and ideal operating table.



Technical Parameters

External Length	2070mm
External Width	550mm
Table Height Min.	700mm
Table Height Max.	1000mm

Technical Configuration

✓ Bed Frame	1pc	✓ Arm Holder with Cushion	2pcs
✓ Mattress	1set	✓ Leg Support with Cushion	2pcs
✓ Anesthesia Screen Frame	1pc	✓ Foot Support with Cushion	2pcs
✓ Shoulder Holder with Cushion	2pcs	✓ Kidney Bridge Handle	1pc

A100-4A

Electric Operating Table
SKL-A

A100-4A is an electric operating table with super load-bearing and extra-long sliding distance, which can be applicable for C-arm. The built-in large-capacity battery can continue to support all operations in the event of a power failure. And the operating table adopts a dual control system to ensure the safety of the operation.



Technical Parameters

External Length	2100mm
External Width	550mm
Table Height Min.	700mm
Table Height Max.	1000mm

Technical Configuration

Bed Frame	1pc	Arm Holder with Cushion	2pcs	Kidney Bridge Handle	1pc
Mattress	1set	Leg Support with Cushion	2pcs		
Anesthesia Screen Frame	1pc	Foot Support with Cushion	2pcs		
Shoulder Holder with Cushion	2pcs	Hand Controller	1pc		

A2000

Electric Operating Table
SKL-A

A2000 is a basic electric operating table with a very high cost effective. It can fulfill the requirement of general surgeries. Its built-in battery can continue to support all operations in the event of a power failure.



Technical Parameters

External Length	2100mm
External Width	480mm
Table Height Min.	680mm
Table Height Max.	930mm

Technical Configuration

Bed Frame	1pc	Arm Holder	2pcs	Sliding Block	10pcs
Mattress	1set	Leg Support	2pcs	Spanner	4pcs
Anesthesia Screen Frame	1pc	Foot Support	2pcs		
Shoulder Holder	2pcs	Hand Controller	1pc		

A205
Manual Operating Table
SKF-C



Technical Parameters

External Length	2000mm
External Width	480mm
Table Height Min.	720mm
Table Height Max.	960mm

Technical Configuration

✓ Arm support	2 pcs	✓ Anaesthetic Screen	1 pc
✓ Leg support	2 pcs		
✓ Body Support	2 pcs		
✓ Shoulder Support	2 pcs		

A3001
Manual Operating Table
SKF-C



Technical Parameters

External Length	2000mm
External Width	480mm
Table Height Min.	790mm
Table Height Max.	1040mm

Technical Configuration

✓ Anesthesia Screen Frame	1 pc	✓ Leg Support with Cushion	2 pcs
✓ Mattress	1 set	✓ Shoulder Support with Cushion	2 pcs
✓ Kidney Bridge Handle	1 pc	✓ Body Support with Cushion	2 pcs
✓ Arm Support with Cushion	2 pcs		



A99-7
Electric Obstetric Bed
SKL-A

250mm
Height
Adjustment

55°/18°
Back Plate
up/down

35°/5°
Seat Plate
up/down

90°
Leg Plate
Swing Out

ELECTRIC OBSTETRIC BED

Frame

Base cover made of high quality ABS. It can be sterilized at the temperature above 100 °C. The tensile strength up to 30 MPa.

Base structure is made of high-quality seamless steel, robot welded with high strength.

Battery

Battery could support about 80-120 complete operation cycles. Working time could be last for one week according to useage and help to keep the operating table running in the event of power failure.



Bed Width
Seamless foam mattress, with waterproof PU cover. Easy cleaning and disinfection.

Arm Plate
Flexible arm plate adjustment.

Hand Controller
Electrical braking system. Controlled by hand controller.

Technical Parameters

External Length	1980mm
External Width	720mm
Table Height Min.	630mm
Table Height Max.	880mm

Technical Configuration

✓ Bed Frame	1pc	✓ Hand Bars	2pcs	✓ Caster	1pc
✓ Foot Pedal	1pc	✓ Leg Support with Cushion	2pcs		
✓ Hand Controller	1pc	✓ Arm Support with Cushion	2pcs		
✓ Anesthesia Screen	1pc	✓ Stainless Steel Basin	1pc		



A99-13
Electric Obstetric Bed
SKL-A

250mm
Height
Adjustment

60°
Back Plate
up

22°
Reverse
Trendelenburg

22°
Trendelenburg

90°
Leg Plate
Swing Out

16°
Seat Plate
up



Material
Base structure is made of high-quality 304 stainless steel, robot welded with high strength and modern appearance.



Mattress
Polyurethane foam mattress, waterproof and seamless, easy cleaning and disinfection.



LED Examination Light
Standard accessory.
More convenient to do gynecological examination.



Arm Plate
Flexible arm plate adjustment.



Battery Could Support
Battery could support about 80-120 complete operation cycles. Working time could be last for one week according to useage and help to keep the operating table running in the event of power failure.



Auxiliary plate
Auxiliary plate is easy for rotating and dismantling.
Rotating angle: 90°

Technical Parameters

External Length	2100mm
External Width	600mm
Table Height Min.	750mm
Table Height Max.	1000mm

Technical Configuration

✓ Armrest	2pcs	✓ SST Rinsing Basin	1pc
✓ Knee Cush	2pcs	✓ Pedal Controller	1pc
✓ Pull Handle	2pcs	✓ LED Examination Light	1pc
✓ Anesthesia Frame	1pc		

A99-5

Electric Obstetric Bed
SKL-A



Technical Parameters

External Length	1880mm
External Width	660mm
Table Height Min.	680mm
Table Height Max.	940mm

Technical Configuration

✓ Arm Support	2pcs	✓ Tightener	2pcs	✓ Slider	4pcs
✓ Handle	2pcs	✓ Hand Controller	1pc		
✓ Leg Support	2pcs	✓ Foot Switch	1pc		
✓ Stirrup	2pcs	✓ Stainless Steel Basin	1pc		

A105

Manual Obstetric Table
SKF-C



Technical Parameters

External Length	1900mm
External Width	600mm
Table Height Min.	700mm
Table Height Max.	900mm

Technical Configuration

✓ Bed Frame	1pc	✓ Leg Support with Cushion	2pcs
✓ Mattress	1pc	✓ Arm Support with Cushion	2pcs
✓ Anesthesia Screen Frame	1pc	✓ Stainless Steel Basin	1pc
✓ Hand Bars	2pcs		

A98-4

Electric Obstetric Bed
SKL-A



Technical Parameters

External Length	1920mm
External Width	920mm
Table Height Min.	650mm
Table Height Max.	900mm

Technical Configuration

✓ Mattress	1pc	✓ Slider	2pcs	✓ Castors	4pcs
✓ Handle	2pcs	✓ Tightener	2pcs	✓ Central Locking Pedal	2pcs
✓ Leg Support	2pcs	✓ Hand Controller	1pc		
✓ ABS Side Rail	2pcs	✓ Stainless Steel Basin	1pc		

A99-8

Electric Gynecological Exam Couch
SKL-B



Technical Parameters

External Length	1400mm
External Width	740mm
Table Height Min.	610mm
Table Height Max.	1100mm

Technical Configuration

✓ Hand Bars	2pcs
✓ Leg Support with Cushion	2pcs
✓ Foot Pedal	1pc
✓ Stainless Steel Basin	1pc

A99-6

Electric Gynecological Exam Couch
SKL-B



Technical Parameters

External Length	1200mm
External Width	500mm
Table Height Min.	500mm
Table Height Max.	1000mm

Technical Configuration

✓ Simple Side Rail	1set
✓ Leg Holder	1set
✓ 5" Castors 4	4pcs
✓ Foot Pedal	1pc

LUX/series

LED
OPERATION
LAMPLUX
LED
operation lamp

Optical System

The bulbs are imported from Germany OSRAM brand, using a special optical system, 360° uniformly illuminated on the observed object, no shadows, true reproduction of the organization color, and achieve constant color temperature under any light intensity conditions. The material is hard and no mercury pollution. Also, the emitted light does not contain infrared and ultraviolet radiation pollution.

Color Rendering Index


The color rendering index reaches 96, which truly reflects the color of human tissues. It is suitable for various surgical scenes and greatly reduces the visual fatigue of medical staff due to long-term operations.



Color Temperature&illumination

Using a new type of LED cold light source, the illuminance can reach 140000 Lux, which can realize 16-level brightness adjustment; the color temperature is in the range of 3500-5000K, and four-level color temperature adjustment can be achieved. While meeting the high illuminance, adjust the lighting parameters according to the needs of different surgeons, so that the surgeon's perception of light is soft and no glint.

OSRAM

CRI
≥96Lux
140000Kelvin
3500-5000

Lens

Each lens adopts a 40mm super large lens, and the lens arrangement is arranged according to the biology of insect compound eyes, and the mutual distance is less than 1mm, so that the shadowless effect is better.

Super Large Lens
40mm

Temperature Rise (Surgical Field)

LED cold light source has no infrared radiation, coupled with our nano-coated radiator, created our excellent heat dissipation effect. It solves the problem of the doctor's head temperature being too high during the operation, and creates a more comfortable operation environment for the doctor; it also effectively avoids the loss of water and infection of the surgical wound, which is conducive to the rapid recovery of the patient's wound after the operation.

Temperature Rise (Surgical Field)
≤1℃

Lighting Mode

With a variety of lighting modes, one-key switching, adapt to different clinical surgery scenes.

Lighting Mode
4

LUXseries
LED Operation Lamp

The LED surgical shadowless lamp is composed of multiple lamp bulbs. The bulb has a longer life than the traditional shadowless lamp. It is round or petal-shaped. It is fixed on the balance arm suspension system. It is stable in positioning and can be moved vertically or cyclically. It can meet different heights during surgery.



The design of the lamp holder focuses on reducing the turbulence. The streamlined lampshade of the whole lamp arm is carefully designed according to the principle of aerodynamics, which is fully in line with the standard of a modern new operating room.

The coating is sprayed with world-leading Akzo Nobel powder coating and anti-glare treatment, which meets the demanding requirements of surgical lighting.



SK-LLD0503E



SK-LLD05E



SK-LLD03E



SK-LLY03E



Optional battery

		SK-LLD0503E	SK-LLD05E	SK-LLD03E	SK-LLY03E
Diameter of lighthouse	mm	860/700	860	700	700
Illumination	Lux	140000/120000	140000	120000	120000
Color Temperature	Kelvin	3500-5000	3500-5000	3500-5000	3500-5000
Color Rendering Index(CRI)	Ra	≥96	≥96	≥96	≥96
Temperature rise(Surgeon head)	°C	≤3	≤3	≤2	≤2
Temperature rise(Surgical field)	°C	≤1	≤1	≤1	≤1
Effective Focusing Depth	cm	100	100	100	100
Lighting Spot Size	cm	22-31	22-31	22-31	22-31
L1+L2 depth of illumination	cm	100	100	100	100
Average Bulb Life	hours	≥50000	≥50000	≥50000	≥50000

Illumination/Lighting Spot Size: at 1m

LUX/series

LED
OPERATION
LAMPLUX
LED
Operation Lamp

Superior Laminar Flow Applications

The streamlined hollow lamp body design reduces the area of the lamp cap that blocks the flow of purified air, thus reducing the air turbulence under the lamp that affects purification.



Green Light Mode

Applicable for surgery requiring high color rendering index, used for cardiothoracic surgery, neurosurgery, etc. Green light can reduce the flare caused by strong light on the eyes.



Special Overhaul Unit

Photoelectric separation design, is for quick repair, which keeps the circuit board away from the heat source. What's more, it increases the service life, and improve the convenience at the same time. The control box can also be easily disassembled for maintenance. The modular lamp body is made of high-quality aluminum alloy precision casting, and the light source can be replaced independently.

One-click Low Light (Endoscope) 3000lux

The operating light is equipped with a one-button low-light (endoscope) illumination mode, which allows the surgeon to quickly switch from regular illumination mode. Provides soft lighting to the environment without interfering with the surgeon's field of view and tissue viewing during the surgical procedure.

Low Light Mode

3000 Lux

Bright

16 Level

Lux

160000

Light Mode

4

Kelvin

3500-5000

High-density Light Mixing System

Using the original high-density light mixing technology, it not only provides ultra-high brightness, but also has high light quality, 16-level brightness adjustment, illuminance 40000-160000Lux, 4-level color temperature adjustment 3500K-5000K, 16-level electric spot size adjustment, Another cavity mirror, red light compensation, green light mode (GL) shadow compensation.

LUXseries
LED Operation Lamp

The LED surgical shadowless lamp is composed of multiple lamp bulbs. The bulb has a longer life than the traditional shadowless lamp. It is round or petal-shaped. It is fixed on the balance arm suspension system. It is stable in positioning and can be moved vertically or cyclically. It can meet different heights during surgery.



The design of the lamp holder focuses on reducing the turbulence. The streamlined lampshade of the whole lamp arm is carefully designed according to the principle of aerodynamics, which is fully in line with the standard of a modern new operating room.

The coating is sprayed with world-leading Akzo Nobel powder coating and anti-glare treatment, which meets the demanding requirements of surgical lighting.



SK-LLD0404D



SK-LLD0303D



SK-LLD04D



SK-LLD03D



		SK-LLD0404D	SK-LLD0303D	SK-LLD04D	SK-LLD03D
Diameter of lighthouse	mm	700/700	615/615	700	615
Illumination	Lux	160000/160000	140000/140000	160000	140000
Color Temperature	Kelvin	3500-5000	3500-5000	3500-5000	3500-5000
Color Rendering Index(CRI)	Ra	≥96	≥96	≥96	≥96
Temperature rise(Surgeon head)	℃	≤2	≤2	≤2	≤2
Temperature rise(Surgical field)	℃	≤1	≤1	≤1	≤1
Effective Focusing Depth	cm	100	100	100	100
Lighting Spot Size	cm	10-25	10-25	10-25	10-25
L1+L2 depth of illumination	cm	100	100	100	100
Average Bulb Life	hours	≥50000	≥50000	≥50000	≥50000

Illumination/Lighting Spot Size: at 1m



IXUS series

LED
OPERATION
LAMP

SK-LLD7272C



SK-LLD7252C



SK-LLD72C



SK-LLD52C



SK-LLY72C



SK-LLY52C



		SK-LLD7272C	SK-LLD7252C	SK-LLD72C	SK-LLD52C	SK-LLY72C	SK-LLY52C
Diameter of lighthead	mm	720/720	720/520	720	520	720	520
Illumination	Lux	160000/160000	160000/160000	160000	160000	160000	160000
Color Temperature	Kelvin	3000-6700	3000-6700	3000-6700	3000-6700	3000-6700	3000-6700
Color Rendering Index(CRI)	Ra	≥98	≥98	≥98	≥98	≥98	≥98
Temperature rise(Surgeon head)	℃	≤1	≤1	≤1	≤1	≤1	≤1
Temperature rise(Surgical field)	℃	≤2	≤2	≤2	≤2	≤2	≤2
Effective Focusing Depth	cm	130	130	130	110	130	110
Lighting Spot Size	cm	10-35	10-35	10-35	10-35	10-35	10-35
L1+L2 depth of illumination	cm	100	100	100	100	100	100
Average Bulb Life	hours	≥60000	≥60000	≥60000	≥60000	≥60000	≥60000

Illumination/Lighting Spot Size: at 1m

WILLYS/series

LED
OPERATION
LAMP

SK-LLD7070A



SK-LLD7050A



SK-LLD70A



SK-LLD50A



SK-LLY70A



SK-LLY50A



		SK-LLD7070A	SK-LLD7050A	SK-LLD70A	SK-LLD50A	SK-LLY70A	SK-LLY50A
Diameter of lighthead	mm	700/700	700/500	700	500	700	500
Illumination	Lux	180000/180000	180000/160000	180000	160000	180000	160000
Color Temperature	Kelvin	3700-5000	3700-5000	3700-5000	3700-5000	3700-5000	3700-5000
Color Rendering Index(CRI)	Ra	85-95	85-95	85-95	85-95	85-95	85-95
Temperature rise(Surgeon head)	℃	≤1	≤1	≤1	≤1	≤1	≤1
Temperature rise(Surgical field)	℃	≤1	≤1	≤1	≤1	≤1	≤1
Effective Focusing Depth	cm	60-130	60-130	60-130	60-130	60-130	60-130
Lighting Spot Size	cm	16-28	16-28	16-28	16-28	16-28	16-28
L1+L2 depth of illumination	cm	130	130	130	130	130	130
Average Bulb Life	hours	≥80000	≥80000	≥80000	≥80000	≥80000	≥80000

Illumination/Lighting Spot Size: at 1m

IXUS/series

LED
OPERATION
LAMP

SK-LLD0505C



SK-LLD05C



SK-LLY05C




SK-LLY03C



		SK-LLD0505C	SK-LLD05C	SK-LLY05C	SK-LLY03C
Diameter of lighthead	mm	650/650	650	650	650
Illumination	Lux	180000/180000	180000	180000	160000
Color Temperature	Kelvin	4000-6000	4000-6000	4000-6000	4000-6000
Color Rendering Index(CRI)	Ra	85-95	85-95	85-95	85-95
Temperature rise(Surgeon head)	°C	≤1	≤1	≤1	≤1
Temperature rise(Surgical field)	°C	≤1	≤1	≤1	≤1
Effective Focusing Depth	cm	60-130	60-130	60-130	60-130
Lighting Spot Size	cm	16-28	16-28	16-28	16-28
L1+L2 depth of illumination	cm	≥130	≥130	≥130	≥130
Average Bulb Life	hours	≥80000	≥80000	≥80000	≥80000

Illumination/Lighting Spot Size: at 1m

WILLYS / series
LED Operation Lamp



		SK-LLD0404B	SK-LLD04B	SK-LLY04B
Diameter of lighthouse	mm	675/675	675	675
Illumination	Lux	160000/160000	160000	160000
Color Temperature	Kelvin	3000-6700	3000-6700	3000-6700
Color Rendering Index(CRI)	Ra	≥97	≥97	≥97
Temperature rise(Surgeon head)	℃	≤1	≤1	≤1
Temperature rise(Surgical field)	℃	≤2	≤2	≤2
Effective Focusing Depth	cm	130	130	130
Lighting Spot Size	cm	10-30	10-30	10-30
L1+L2 depth of illumination	cm	≥120	≥120	≥120
Average Bulb Life	hours	≥60000	≥60000	≥60000

Illumination/Lighting Spot Size: at 1m

VOLTH / series
Integral Reflection Operation Lamp



		SK-LZD7070A	SK-LZD7050A	SK-LZD70A	SK-LZD50A	SK-LZY70A	SK-LZY50A
Diameter of lighthouse	mm	700/700	700/500	700	500	700	500
Illumination	Lux	180000/180000	180000/160000	180000	160000	180000	160000
Color Temperature	Kelvin	4500±250	4500±250	4500±250	4500±250	4500±250	4500±250
Color Rendering Index(CRI)	Ra	≥85	≥85	≥85	≥85	≥85	≥85
Temperature rise(Surgeon head)	℃	≤1	≤1	≤1	≤1	≤1	≤1
Temperature rise(Surgical field)	℃	≤1	≤1	≤1	≤1	≤1	≤1
Effective Focusing Depth	cm	60-120	60-120	60-120	60-120	60-120	60-120
Lighting Spot Size	cm	16-28	16-28	16-28	16-28	16-28	16-28
L1+L2 depth of illumination	cm	≥120	≥120	≥120	≥120	≥120	≥120
Average Bulb Life	hours	≥1500	≥1500	≥1500	≥1500	≥1500	≥1500

Illumination/Lighting Spot Size: at 1m

SAIKANG


MEDICAL BRIDGE CEILING PENDANT



SK-PQ003

Combined Medical Bridge Ceiling Pendant

The combined medical bridge ceiling pendant is suitable for compact rooms, two or more modules of medical bridge ceiling pendant can be connected from head to end, improving the space utilization of ICU room, saving treatment space as much as possible without all configurations changed.




SK-PQ001

Medical Bridge Ceiling Pendant (separate wet and dry areas)

The ICU medical bridge ceiling pendant (separate wet and dry areas) is a medical rescue auxiliary equipment.

Features:


1. There are multi-layer instrument platforms in wet (gas) and dry (electric) areas;
2. The height of drawers can be adjusted;
3. With retractable IV pole, with four hocks;
4. The function columns and cabinets built in rotation limit device to prevent the drift;



SK-PQ002

Medical Bridge Ceiling Pendant (Integration of dry and wet areas)

The medical bridge ceiling pendant (integration of dry and wet areas) is based on the separation of dry and wet, the two areas are combined into one, suitable for NICU and small nursing areas.



SK-PQ004

Multi-arm Medical Ceiling Pendant

The multi-arm combined medical bridge ceiling pendant is suitable for compact ICU rooms, the combination of flexible arms and function columns can make full use of the space of the ICU rooms.



SK-PQ005

Multi-arm Medical Ceiling Pendant

The multi-arm combined medical bridge ceiling pendant is suitable for compact ICU rooms, the combination of flexible arms and function columns can make full use of the space of the ICU rooms.



SK-PT001

Single Arm Medical Bridge Ceiling Pendant

The single arm medical bridge ceiling pendant, meeting the daily load-bearing requirements and using strong material, made of imported antibacterial plastic powder, ensuring the cleanliness of the operating environment. Its rich accessories matching various systems, greatly extending the clinical practicability.



SK-PT002

Single Arm Medical Endoscopic Bridge Ceiling Pendant

The single arm medical endoscopic bridge ceiling pendant is endoscopy surgery auxiliary equipment, meeting the daily load-bearing requirements and using strong material, made of imported antibacterial plastic powder, ensuring the cleanliness of the operating environment. Its rich accessories matching various systems, greatly extending the clinical practicability.






SK-PT003

Double Arms Medical Bridge Ceiling Pendant

The double arms medical bridge ceiling pendant is more flexible in moving and positioning, compared with the single one, longer treatment coverage radius, meeting the daily load-bearing requirements and using strong material, made of imported antibacterial plastic powder, ensuring the cleanliness of the operating environment. Its rich accessories matching various systems, greatly extending the clinical practicability.



SK-PT004

Mobile Medical Bridge Ceiling Pendant

The mobile medical bridge ceiling pendant can be connected with external power supply and medical air source, bearing 200kg, prevent the shortages of gas electricity supply when adding nursing machines.



SK-PT005

Anesthesia Medical Bridge Ceiling Pendant

The anesthesia medical bridge ceiling pendant is equipped with hanging holder for anesthesia machine, bearing 250 kg of anesthesia machine, with international AGSS jet exhaust gas, suitable for anesthesia surgery.



SK-PT006

Supply Column


Working table supports can be placed on both sides of the supply column, one-stop operation for medical staff, suitable for supply rooms.



SK-PT008

Medical Bridge Ceiling Pendant For Imaging Operating Room

The medical medical bridge ceiling pendant for imaging operating room is a basic equipment specially designed for endoscopic treatment. Used for carrying instrument of endoscope, monitor and so on. Providing air pump, power supply, video and communication interfaces for equipment, with cable and pipeline, making the operating environment orderly.



ACCESSORIES

Operating Table Accessories

Floor Type Multi-functional Orthopedic Frame



Connecting Multi-functional Orthopedic Frame



Carbon Fiber Multi-functional Orthopedic Frame



Multi-purpose Head Rack



Shoulder Traction Frame



Mechanical Power IV Pole



Stirrup Pneumatic Power Leg Rack



Bow Spine Operating Supporter



Operating Table Accessories Cart



Knee surgery rack



Double table top



Prone body supporting frame



Shoulder surgery rack



Arm microsurgical rack



Bed surface widening module



External kidney bridge



Dirt basin



Body safe strap



ACCESSORIES

Gel Mattress Position Solution

Supine Position



Prone Position









Prone Position




Cut Stone Position



Medical Bridge Ceiling Pendant Accessories

Disposable Oxygen Inhaler II	Disposable Oxygen Inhaler A	Disposable Oxygen Inhaler (Simple Type)
		
External Hanging Oxygen Inhaler	External Hanging Waste Liquid Collection Device	Oxygen Inhaler (All Metal)
		

Waste Liquid Collection Device	Negative Pressure Regulating Valve	Standard Medical Terminals
		
Interface Medical Terminal	AGSS Jet Exhaust Device	Terminal Plug
		



SERVICE NETWORK

Stable product quality and reliable supply reputation have won the favor of well-known enterprises at home and abroad. The company's products are all over Europe, America, Southeast Asia and other regions, facing the world. With the modern information-based sales network, the company strides into the fast lane of professional and large-scale development. The company looks forward to working with you to create a better future.



Saikang has delivered more than
1000000
sets of products to the world




7×24
our service hotline,
always responsive to
customer

MORTUARY EQUIPMENT



SK-DT004

AUTOPSY TABLE



Technical parameters


▪ External size (LxWxH)	2500x800x850 mm
▪ Safe working load	150 kg

Technical configuration

<input checked="" type="checkbox"/> Sink	1pc
<input checked="" type="checkbox"/> Socket safety	1pc
<input checked="" type="checkbox"/> Water-tap	2pcs
<input checked="" type="checkbox"/> Water spray gun	1pc

Optional

<input type="checkbox"/> Height adjustable table top	
<input type="checkbox"/> Three punching plates	
<input type="checkbox"/> Crush machine	
<input type="checkbox"/> Pillow	
<input type="checkbox"/> Intrusment tray	



Material

Baosteel from Fortune Global 500. Panasonic Robotic ensure 360° full smooth welding.

**Table top**

304 stainless steel table top, with two drainage through to ensure the water fluids to the drain easily. The table top with the rounded corners, easy to clean.

**Sink**

With a hair isolator in the sink to prevent the drain pipe from being blocked.

**Washing machine**


Equipped a manual and a sensible water-tap, with a heat and cold regulator, effectively prevent cross-infection of virus. (need to connect a water heater)

**Water spray gun**

With a water spray gun, adjustable water pressure, easy to flush the entire dissection table.

**Ventilation system**

Both the inner side of the table are inclined downward to exhaust the harmful gas during the dissection process. (need to connect the external ventilation system)

**Socket safety**

220V.50HZ waterproof socket, with a leakage protection device.

**Crush**

Optional accessory: with the crush machine to prevent the drain pipe from being blocked.

**Optional function**

With the height adjustable table top:750-950mm.

**Optional accessory**


Three punching plates, easy to disassemble and clean.

**Optional accessory**

With the pillow, more stable to fix the head.

**Optional accessory**

Stainless steel Intrusment tray.



SK-DT003

AUTOPSY TABLE

Technical parameters

▪ External size (LxWxH)	2400x700x850 mm
▪ Safe working load	150 kg

Technical configuration

<input checked="" type="checkbox"/> Drain pipe	1pc
<input checked="" type="checkbox"/> Castors	4pcs



Material

High quality 304 stainless steel material with mirror surface polishing. Panasonic Robotic ensure 360° full smooth welding.



Table top

Environmental plastic can be sterilized in 100°C, tensile upto 30MPa.



Disassembly structure

Disassembly structure, save volume, easy to be stored.



Castors

Four silent casters of 3 inch diameter, two with brakes, easy to install, stable and smooth.



Sink

Size of the sink:
(L*W*H)450*320*180mm
With a hair isolator in the sink to prevent the drain pipe from being blocked.



Drain pipe

With the adjustable drain pipe, easy to install.



SK-DT005

AUTOPSY TABLE



Technical parameters

▪ External size (LxWxH)	2050x820x(720-1040) mm
▪ Safe working load	250 kg
▪ Reverse trendelenburg	0-15°
▪ Trendelenburg	0-15°

Technical configuration

<input checked="" type="checkbox"/> USA Hydraulic Pump	2pcs
<input checked="" type="checkbox"/> Hydraulic Foot Pedals	3pcs
<input checked="" type="checkbox"/> Central Locking Pedal	2pcs
<input checked="" type="checkbox"/> Castors	4pcs



Metal Material

Baosteel from Fortune Global 500. Panasonic Robotic ensure 360° full smooth wedding.



Plastic

Environmental plastic can be sterilized in 100°C, tensile upto 30MPa.



Table Top

Stainless steel 304 table top, deep trough around all 4 sides assists the easy flow of body fluids to the drain at the center of the foot end.



Hydraulic Pumps

3 hydraulic foot pedals make raising or lowering either or both ends of the table.



Brake

Dual side mounted foot controls allow for almost effortless operating on either side. Dual locking levers lock all the wheels simultaneously from either end.



Castors

Four double side casters of 8 inch diameter, simultaneously lockable by the step pedal at both sides.



SKB-7C011

HYDRAULIC CORPSE TRANSPORTATION TROLLEY

Technical parameters

▪ External size (LxWxH)	2200x670x400-1700 mm
▪ Body tray size (LxW)	1950x610 mm
▪ Height adjustment	400-1700 mm
▪ Safe working load	400 kg

Technical configuration

<input checked="" type="checkbox"/> Body tray	1 pc
<input checked="" type="checkbox"/> Castors	4 pcs
<input checked="" type="checkbox"/> Hydraulic cylinder	1 pc
<input checked="" type="checkbox"/> Pedal	1 pc

Material

High quality 304 stainless steel material with mirror surface polishing. Panasonic Robotic ensure 360° full smooth welding. Single-fork structure, durable and stable.




Body tray

Top of the elevator equipped with a roller tray, lockable and easy to move.




Size

The size of the body tray could be customized.



Height adjustment

Height adjustment is controlled by the pedal hydraulic cylinder. The pedal surface with anti-skid treatment.




Castors

With two directional castors and two universal castors, stable and smooth.



SKB-7C011-1

HYDRAULIC CORPSE TRANSPORTATION TROLLEY



Technical parameters

▪ External size (LxWxH)	2200x670x450-2000 mm
▪ Body tray size (LxW)	1950x610 mm
▪ Height adjustment	450-2000 mm
▪ Safe working load	300 kg

Technical configuration

<input checked="" type="checkbox"/> Body tray	1 pc
<input checked="" type="checkbox"/> Castors	4 pcs
<input checked="" type="checkbox"/> Hydraulic cylinder	1 pc
<input checked="" type="checkbox"/> Pedal	1 pc




Material

High quality 304 stainless steel material with mirror surface polishing. Panasonic Robotic ensure 360° full smooth welding. Double-fork structure, durable and stable.



Height adjustment

Height adjustment is controlled by the pedal hydraulic cylinder. The pedal surface with anti-skid treatment.



Body tray

Top of the elevator equipped with a roller tray, lockable and easy to move.



Castors

With two directional castors and two universal castors, stable and smooth.




Size

The size of the body tray could be customized.



SKB-7C012

ELECTRIC HYDRAULIC CORPSE TRANSPORTATION TROLLEY



Technical parameters

▪ External size (LxWxH)	2200x670x400-1700 mm
▪ Body tray size (LxW)	1950x610 mm
▪ Height adjustment	400-1700 mm
▪ Safe working load	400 kg

Technical configuration

<input checked="" type="checkbox"/> Body tray	1 pc
<input checked="" type="checkbox"/> Castors	4 pcs
<input checked="" type="checkbox"/> Electric hydraulic cylinder	1 pc
<input checked="" type="checkbox"/> Back up battery	1 pc



Material

High quality 304 stainless steel material with mirror surface polishing. Panasonic Robotic ensure 360° full smooth welding. Single-fork structure, durable and stable.



Height adjustment

Height adjustment by electric hydraulic cylinder with the battery back up. The back up battery can load 50kg body and work for 100 times.



Body tray

Top of the elevator equipped with a roller tray, lockable and easy to move.



Castors

With two directional casters and two universal casters, stable and smooth.




Size

The size of the body tray could be customized.



SKB-7C012-1

ELECTRIC HYDRAULIC CORPSE TRANSPORTATION TROLLEY



Technical parameters

▪ External size (LxWxH)	2200x670x450-2000 mm
▪ Body tray size (LxW)	1950x610 mm
▪ Height adjustment	450-2000 mm
▪ Safe working load	300 kg

Technical configuration

<input checked="" type="checkbox"/> Body tray	1 pc
<input checked="" type="checkbox"/> Castors	4 pcs
<input checked="" type="checkbox"/> Electric hydraulic cylinder	1 pc
<input checked="" type="checkbox"/> Back up battery	1 pc



Material

High quality 304 stainless steel material with mirror surface polishing. Panasonic Robotic ensure 360° full smooth welding. Double-fork structure, durable and stable.



Height adjustment

Height adjustment by electric hydraulic cylinder with the battery back up. The back up battery can load 50kg body and work for 100 times.



Body tray

Top of the elevator equipped with a roller tray, lockable and easy to move.



Castors

With two directional castors and two universal casters, stable and smooth.



Size

The size of the body tray could be customized.



SKB-7A001

MORTUARY REFRIGERATOR

Technical parameters

▪ Rear compressor size (LxWxH)	2280x800x610 mm
▪ Upper compressor size (LxWxH)	2060x800x900 mm
▪ Opening size (LxW)	630x440 mm
▪ Pallet size (LxW)	1950x580 mm
▪ Rated temperature	5~-15°C
▪ Environment temperature	10-38°C
▪ Power supply	AC220W, 50HZ
▪ Refrigerant	R134A

Technical configuration

<input checked="" type="checkbox"/> Door sealing	1 pc
<input checked="" type="checkbox"/> Compressor	1 pc
<input checked="" type="checkbox"/> Body tray	1 pc
<input checked="" type="checkbox"/> LED light	1 pc
<input checked="" type="checkbox"/> Control panel	1 pc
<input checked="" type="checkbox"/> Condenser	1 pc



Material

High quality stainless steel with the Mirror surface polishing.
Panasonic Robotic ensure 360° full smooth welding surface.




Refrigeration

Quick refrigeration, stably keep the temperature within 5°C to -15°C.



Door sealing

One built heavy cabinet door, thickness:50mm.
High quality and undeformable during normal using.
Replaceable door seal close to the surface to make sure stable the temperature.



Refrigeration system

High-strength compressor with good refrigeration and energy saving performance.
Each room has an independent compressor to make sure every room will not influence each other.



Body tray

High quality stainless steel body tray with high stability and load capacity.



Light

Aluminum alloy LED light with longer service life and more stable brightness.



Temperature control

High precision digital display of control panel.
With the sound and light alarm of over temperature.



Condenser

Closely distributed copper tubes, with the strong refrigeration effect.



SKB-7A002

MORTUARY REFRIGERATOR

Technical parameters

▪ Rear compressor size (LxWxH)	2280x800x1300 mm
▪ Upper compressor size (LxWxH)	2060x800x1475 mm
▪ Opening size (LxW)	630x440 mm
▪ Pallet size (LxW)	1950x580 mm
▪ Rated temperature	5~-15°C
▪ Environment temperature	10-38°C
▪ Power supply	AC220W, 50HZ
▪ Refrigerant	R134A

Technical configuration

<input checked="" type="checkbox"/> Door sealing	2 pcs
<input checked="" type="checkbox"/> Compressor	2 pcs
<input checked="" type="checkbox"/> Body tray	2 pcs
<input checked="" type="checkbox"/> LED light	2 pcs
<input checked="" type="checkbox"/> Control panel	2 pcs
<input checked="" type="checkbox"/> Condenser	2 pcs



Material

High quality stainless steel with the Mirror surface polishing.
Panasonic Robotic ensure 360° full smooth welding surface.




Refrigeration

Quick refrigeration, stably keep the temperature within 5°C to -15°C.



Door sealing

One built heavy cabinet door, thickness:50mm.
High quality and undeformable during normal using.
Replaceable door seal close to the surface to make sure stable the temperature.



Refrigeration system

High-strength compressor with good refrigeration and energy saving performance.
Each room has an independent compressor to make sure every room will not influence each other.



Body tray

High quality stainless steel body tray with high stability and load capacity.



Light

Aluminum alloy LED light with longer service life and more stable brightness.



Temperature Control

High precision digital display of control panel.
With the sound and light alarm of over temperature.



Condenser

Closely distributed copper tubes, with the strong refrigeration effect.



SKB-7A003

MORTUARY REFRIGERATOR

Technical parameters

▪ Rear compressor size (LxWxH)	2280x800x1785 mm
▪ Upper compressor size (LxWxH)	2060x800x1975 mm
▪ Opening size (LxW)	630x440 mm
▪ Pallet size (LxW)	1950x580 mm
▪ Rated temperature	5~-15°C
▪ Environment temperature	10-38°C
▪ Power supply	AC220W, 50HZ
▪ Refrigerant	R134A

Technical configuration

<input checked="" type="checkbox"/> Door sealing	3 pcs
<input checked="" type="checkbox"/> Compressor	3 pcs
<input checked="" type="checkbox"/> Body tray	3 pcs
<input checked="" type="checkbox"/> LED light	3 pcs
<input checked="" type="checkbox"/> Control panel	3 pcs
<input checked="" type="checkbox"/> Condenser	3 pcs



Material

High quality stainless steel with the Mirror surface polishing.
Panasonic Robotic ensure 360° full smooth welding surface.




Refrigeration

Quick refrigeration, stably keep the temperature within 5°C to -15°C.



Door sealing

One built heavy cabinet door, thickness:50mm.
High quality and undeformable during normal using.
Replaceable door seal close to the surface to make sure stable the temperature.



Refrigeration system

High-strength compressor with good refrigeration and energy saving performance.
Each room has an independent compressor to make sure every room will not influence each other.



Body tray

High quality stainless steel body tray with high stability and load capacity.



Light

Aluminum alloy LED light with longer service life and more stable brightness.



Temperature Control

High precision digital display of control panel.
With the sound and light alarm of over temperature.



Condenser

Closely distributed copper tubes, with the strong refrigeration effect.



SKB-7A004

MORTUARY REFRIGERATOR

Technical parameters

▪ Rear compressor size (LxWxH)	2280x1580x1300 mm
▪ Upper compressor size (LxWxH)	2060x1580x1475 mm
▪ Opening size (LxW)	630x440 mm
▪ Pallet size (LxW)	1950x580 mm
▪ Rated temperature	5~-15°C
▪ Environment temperature	10-38°C
▪ Power supply	AC220W, 50HZ
▪ Refrigerant	R134A

Technical configuration

<input checked="" type="checkbox"/> Door sealing	4 pcs
<input checked="" type="checkbox"/> Compressor	4 pcs
<input checked="" type="checkbox"/> Body tray	4 pcs
<input checked="" type="checkbox"/> LED light	4 pcs
<input checked="" type="checkbox"/> Control panel	4 pcs
<input checked="" type="checkbox"/> Condenser	4 pcs



Material

High quality stainless steel with the Mirror surface polishing.
Panasonic Robotic ensure 360° full smooth welding surface.




Refrigeration

Quick refrigeration, stably keep the temperature within 5°C to -15°C.



Door sealing

One built heavy cabinet door, thickness:50mm.
High quality and undeformable during normal using.
Replaceable door seal close to the surface to make sure stable the temperature.



Refrigeration system

High-strength compressor with good refrigeration and energy saving performance.
Each room has an independent compressor to make sure every room will not influence each other.



Body tray

High quality stainless steel body tray with high stability and load capacity.



Light

Aluminum alloy LED light with longer service life and more stable brightness.



Temperature Control

High precision digital display of control panel.
With the sound and light alarm of over temperature.



Condenser

Closely distributed copper tubes, with the strong refrigeration effect.



SKB-7A006

MORTUARY REFRIGERATOR

Technical parameters

▪ Rear compressor size (LxWxH)	2280x1580x1785 mm
▪ Upper compressor size (LxWxH)	2060x1580x1975 mm
▪ Opening size (LxW)	630x440 mm
▪ Pallet size (LxW)	1950x580 mm
▪ Rated temperature	5~-15°C
▪ Environment temperature	10-38°C
▪ Power supply	AC220W, 50HZ
▪ Refrigerant	R134A

Technical configuration

<input checked="" type="checkbox"/> Door sealing	6 pcs
<input checked="" type="checkbox"/> Compressor	6 pcs
<input checked="" type="checkbox"/> Body tray	6 pcs
<input checked="" type="checkbox"/> LED light	6 pcs
<input checked="" type="checkbox"/> Control panel	6 pcs
<input checked="" type="checkbox"/> Condenser	6 pcs



Material

High quality stainless steel with the Mirror surface polishing.
Panasonic Robotic ensure 360° full smooth welding surface.



Refrigeration

Quick refrigeration, stably keep the temperature within 5°C to -15°C.



Door sealing

One built heavy cabinet door, thickness:50mm.
High quality and undeformable during normal using.
Replaceable door seal close to the surface to make sure stable the temperature.



Refrigeration system

High-strength compressor with good refrigeration and energy saving performance.
Each room has an independent compressor to make sure every room will not influence each other.



Body tray

High quality stainless steel body tray with high stability and load capacity.



Light

Aluminum alloy LED light with longer service life and more stable brightness.



Temperature Control

High precision digital display of control panel.
With the sound and light alarm of over temperature.



Condenser

Closely distributed copper tubes, with the strong refrigeration effect.



SKB-7A008

MORTUARY REFRIGERATOR

Technical parameters

▪ Rear compressor size (LxWxH)	2400x1580x2110 mm
▪ Opening size (LxW)	630x440 mm
▪ Pallet size (LxW)	1950x580 mm
▪ Rated temperature	5~-15°C
▪ Environment temperature	10-38°C
▪ Power supply	AC220W, 50HZ
▪ Refrigerant	R134A

Technical configuration

<input checked="" type="checkbox"/> Door sealing	8 pcs
<input checked="" type="checkbox"/> Compressor	8 pcs
<input checked="" type="checkbox"/> Body tray	8 pcs
<input checked="" type="checkbox"/> LED light	8 pcs
<input checked="" type="checkbox"/> Control panel	8 pcs
<input checked="" type="checkbox"/> Condenser	8 pcs



Material

High quality stainless steel with the Mirror surface polishing.
Panasonic Robotic ensure 360° full smooth welding surface.




Refrigeration

Quick refrigeration, stably keep the temperature within 5°C to -15°C.



Door sealing

One built heavy cabinet door, thickness:50mm.
High quality and undeformable during normal using.
Replaceable door seal close to the surface to make sure stable the temperature.



Refrigeration system

High-strength compressor with good refrigeration and energy saving performance.
Each room has an independent compressor to make sure every room will not influence each other.



Body tray

High quality stainless steel body tray with high stability and load capacity.



Light

Aluminum alloy LED light with longer service life and more stable brightness.



Temperature Control

High precision digital display of control panel.
With the sound and light alarm of over temperature.



Condenser

Closely distributed copper tubes, with the strong refrigeration effect.



SKB-7B011

BODY BAG

Technical parameters

▪ Size (LxW)	2200x90 mm
▪ Thickness	0.2 mm
▪ Safe working load	200 kg

Technical configuration

<input checked="" type="checkbox"/> Handle rope	6 pcs
<input checked="" type="checkbox"/> Zipper	2 pcs

Optional

☐ Name card holder



Material

Using the PEVA Environmental friendly material.
The surface with waterproof coating to prevent leaking of body fluid.



Handle rope

With six strong handle ropes, easy to transport.



Zipper

Double side of zippers design, easy to open and close.



Name card holder

Optional to add the name card holder, easy to check the identity.



SKH107

CORPSE TRANSPORTATION TROLLEY

Technical parameters

▪ External size (LxWxH)	2050x780x900 mm
▪ Safe working load	150 kg

Technical configuration

<input checked="" type="checkbox"/> Castors	4pcs
---	------

Optional

<input type="checkbox"/> Table cover	
--------------------------------------	--



Material

High quality 304 stainless steel material with mirror surface polishing. Panasonic Robotic ensure 360° full smooth welding.



Table top

Special design of the hollow table top to prevent the dirt from spilling out.



Castors

Four silent casters of 5 inch diameter, all with brakes, separately lockable by the step pedal, stable and smooth.



Table cover

Optional: with the 304 stainless steel table cover.



SKB-7C010


FEATURES AND HIGHLIGHTS

Technical parameters

▪ External size (LxWxH)	1950x730x680-900 mm
▪ Height adjustment	680-900 mm
▪ Safe working load	300 kg

Technical configuration

<input checked="" type="checkbox"/> Body tray	1 pc
<input checked="" type="checkbox"/> Castors	4 pcs
<input checked="" type="checkbox"/> Privacy protection cover	1 pc
<input checked="" type="checkbox"/> Pedal	1 pc



Material

High quality 304 stainless steel material with mirror surface polishing. Panasonic Robotic ensure 360° full smooth welding. Unique lifting structure, durable and stable.



Function

Height adjustment is controlled by the pedal hydraulic, providing strong support for the functions of height adjustment. The pedals surface are with anti-skid treatment.



Body tray

Tray top with removable structure, height 440mm.



Castors

With four silent casters of 8 inch diameter, locked by the step pedal. Flexible steering, easy to move, can easily drive fast on bumpy roads.



Privacy protection cover

Black privacy protection cover, protect the privacy of the deceased, with zipper, removable, water proof.



SKB-7D001

CORPSE STORAGE SHELF

Technical parameters

▪ Size (LxWxH)	2000x620x1800 mm
▪ Body tray size (LxW)	1950x580 mm
▪ Safe working load	150 kg


Technical configuration

<input checked="" type="checkbox"/> Body tray	4 pcs
<input checked="" type="checkbox"/> Castors	4 pcs



Material

High quality 304 stainless steel with the Mirror surface polishing. Panasonic Robotic ensure 360° full smooth welding.



Function

Disassembly structure, save storage and shipping volume.



Body tray

Each layer equipped a body tray with handles, strong and stable. Each shelf is with rollers, lockable and can be smoothly pulled out.



Castors

Four silent castors of 3 inch diameter, separately lockable by the step pedal.



Customize

The number and size of storage layers can be customized according to the customer's space size and requirements.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service

WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa





Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi,
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SKMortech™
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
 +86 18021231901 ·  24/7
 export@skmortech.com
 www.skmortech.com



Company video



General catalogue